PLASMAFILTER

INSTRUCTIONS FOR USE

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(ENGLISH)

(CHINESE)

(TAIWANESE)

Gambro Dialysatoren GmbH
Holger-Crafoord-Strasse 26
72379, Hechingen, Germany

N50 184; version 006
Date of issue: 2015-05-01

© Copyright Gambro Dialysatoren GmbH

 **GAMBRO**

General information

The information contained in this instruction manual must be read before the use of this device to ensure a safe and effective treatment.

The plasmafilter is delivered sterile with non-pyrogenic fluid compartments and is ready for use after priming and rinsing.

The plasmafilter is for single use only. The quality of the product is guaranteed only when it is prepared and used according to the procedures described in the "Directions for use".

Recommended storage conditions: Below + 30°C (+ 86°F).

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

WARNING ! is used to alert the operator not to take a specific action that can cause potential hazard which, if ignored, can result in an adverse reaction, injury or death.

CAUTION ! is used to alert the operator to take a specific action to protect against possible hazard which, if ignored, could have an adverse effect on the patient or equipment.

NOTE ! is used to remind the operator of normal treatment functions and what is the suitable action in a particular situation.

ADVERSE REACTIONS

Certain types of adverse reactions may occur due to operational factors associated with the treatment. Therefore, proper management of the fluid removal, anticoagulation, blood, ultrafiltrate and infusion solution flow rates as well as monitoring of the overall treatment parameters are essential to avoid side effects which may be associated with plasmafiltration.

Hypersensitivity reactions have been observed. Should such symptoms appear, particularly at onset of the treatment session, it is important to react immediately by discontinuing the session and administering appropriate treatment. Blood in the extracorporeal circuit must not be returned to the patient. Extra care must be taken when treating patients who have exhibited possible hypersensitivity symptoms during previous treatments, or patients with known or suspected hypersensitivity to products sterilized with ethylene oxide, or patients who have a history of being highly sensitive and allergic to a variety of substances. A physician must be consulted to evaluate the risk and prescribe the appropriate precautions if a possible sensitivity is suspected.

The following factors are considered essential to minimize the risk of a hypersensitivity reaction and other side effects:

- The choice of plasmafilter configuration as well as type of membrane and sterilization method/agent.
- Strict adherence to the set-up, priming and rinsing procedures detailed in the manufacturer's "Directions for use".
- Setting and monitoring the treatment operating parameters according to the manufacturer's recommendations specified for each type of plasmafilter and to the patient's needs and tolerance.
- Strict adherence to all WARNINGS and CAUTIONS given by the manufacturer in the "Directions for use".

GENERAL WARNINGS AND CAUTIONS

Do not perform the treatment using operating values outside the range specified for the plasmafilter.

Do not use chemicals or solvents neither inside nor outside the plasmafilter without approval from the manufacturer. Halogenated solvents (perchloroethylene, chloroform, etc.), halogen containing liquids (e.g. iodine tinctures), denatured ethyl alcohol, acetone, phenols or their derivatives must be avoided.

Blood and ultrafiltrate connectors on this plasmafilter are designed to comply with international and national standards. Use only with blood lines and fluid connectors that comply with international and national standards.

The risk of cross contamination between patients may be

prevented by the use of pressure transducer protectors.

Patients receiving extracorporeal treatment such as plasmafiltration can be exposed to the risk of infection, therefore aseptic technique is required throughout the entire set-up and treatment procedures.

Patients on life support equipment should be visually monitored by competent medical personnel since hazardous circumstances can arise which may not activate alarms. Constant monitoring of all relevant parameters must be maintained throughout the treatment session to reduce the risk of adverse reactions and clinical side effects (e.g. clotting, blood loss and air embolism). Adjust all operating parameters to prescribed value range indicated by the physician for the specific patient's needs and tolerance level.

WARNING ! Processing this dialyzer on a reuse machine or manual reuse processes may cause serious damage to the product resulting in blood leaks.

CAUTION ! Federal law (USA) restricts this device to sale by or on the order of a physician.

WARRANTY AND LIMITATION OF LIABILITY

The manufacturer warrants that the plasmafilter has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the plasmafilter becoming apparent before the expiry date.

Expiry date: Refer to information on the unit container label.

The warranty above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, expressed or implied, statutory or otherwise, and there are no warranties of merchantability or fitness or other warranties, which extend beyond those described.

The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the plasmafilter and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the plasmafilter, whether as a result of any defect therein or otherwise.

The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings, directions and instructions in the labelling, damage arising from events after the manufacturer's release of the plasmafilter, failure or omission to inspect the plasmafilter before use in order to insure that the plasmafilter is in proper condition, or any warranty given by independent distributors or dealers.

CLINICAL SIDE EFFECTS

If clotting occurs in the vascular access devices stop the blood pump, clamp the blood lines and disconnect the blood lines from the patient. Unclamp the lines and start the blood pump to allow the blood to circulate through the plasmafilter with some positive pressure (30 - 100 mmHg) in the system, e.g. by means of a clamp. Stop the infusion and filtrate pumps and clamp the filtrate line. Check the clotting time of the circulating blood and administer heparin, if necessary. Decloit the vascular access devices, if necessary, replace the needles. Stop the blood pump, clamp the lines and connect them to the patient. Unclamp the lines and start the blood pump, unclamp the filtrate line and start the infusion and filtrate pumps to restart the treatment.

CAUTION ! It is recommended not to return the blood to the patient if the blood has been recirculated too long, to avoid risk for adverse reactions.

If clotting occurs in the plasmafilter both the plasmafilter and the blood lines must be replaced. Stop the blood pump, clamp the lines and disconnect the patient. Flush the vascular access devices with heparin and physiological solution according to the procedure used in the clinic. Discard the plasmafilter and the blood lines.

WARNING ! Do not return the blood to the patient.

If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the plasmafilter. If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

If an internal blood leakage may occur initiate appropriate action. The operator, upon confirming a blood leak, must stop the treatment session and replace the plasmafilter.

If air enters the extracorporeal blood circuit air embolism is caused and can be hazardous to the patient. To minimize the risk of air embolism constant monitoring of the extracorporeal blood circuit, both visually and with an air detector, is recommended.

Strict adherence to the manufacturer's recommended pre-treatment activities/set-up, set-up and start-up procedures will facilitate the removal and prevent accumulation of air in the plasmafilter **before** the treatment session. If air enters or is identified in the plasmafilter during priming and rinsing and cannot be removed with the use of additional physiological solution, the plasmafilter must be replaced.

WARNING ! Do not initiate or continue the treatment session if air that can cause a hazard to the patient is present in the extracorporeal circuit. Accidental disconnection or air leakage may cause fatal air embolism.

Immediately stop the treatment session if air is identified in the blood circuit. Discontinue the treatment session and do not attempt to return the blood to the patient if the air cannot be safely removed.

If a hypersensitivity reaction appears discontinue the treatment session immediately by stopping the blood pump, clamping the blood lines and disconnecting the patient. Initiate appropriate medical intervention. Discard the plasmafilter and the blood lines.

WARNING ! The blood in the extracorporeal circuit must not be returned to the patient.

N50 008 / 004

Directions for use

CAUTION ! Strict adherence to the procedural order and content, including all operational WARNINGS and CAUTIONS, described in the following procedures by the manufacturer is required to ensure a safe and effective treatment.

CAUTION ! Aseptic technique is required throughout plasmfiltration preparation and treatment to avoid contamination.

PRE-TREATMENT SET UP

CAUTION ! Always, check the plasmafilter carefully and use only if it is mechanically intact and the protection caps over the blood and filtrate ports are securely in place, dry and intact. Do not use the plasmafilter if the unit (or sterile) package is damaged because sterility may be compromised.

Position the plasmafilter vertically in the appropriate holder with the arterial end at the bottom (1).

Remove the protection caps from the blood ports on the plasmafilter and attach the blood line connectors; the arterial line at the bottom of the plasmafilter and the venous line at the top.

Remove the cap with medical grade paper from the filtrate port and attach the filtrate line to the open filtrate port of the plasmafilter (2). Remove the filtrate pump segment from the filtrate pump housing and clamp the filtrate line close to the plasmafilter.

WARNING ! Do not remove the dense screw cap from the second port of the plasmafilter.

Place clamps on the blood lines according to the procedure followed in the unit.

Prepare and connect a physiological solution to the arterial blood line for priming and rinsing the plasmafilter:

- PF 1000 N at least 2000 ml
- PF 2000 N at least 3000 ml

Place the free end of the venous line in an empty container.

Connect the heparin line to a heparin container.


START UP, PRIMING AND RINSING

A minimum of priming solution is required to rinse and prime the extracorporeal circuit.


- PF 1000 N = 2000 ml
- PF 2000 N = 3000 ml

The last 1000 ml should contain 5000 IU heparin.

Remove any obstructing clamps and start the infusion pump to rinse and prime the infusion circuit.

Remove any obstructing clamps and start the blood pump at 50 ml/min (3) to rinse and prime the blood circuit. Increase the flow rate to 100 ml/min. When the solution reaches the venous line, clamp the filtrate side (4).

After 500 ml, open the filtrate line and clamp the venous line (5). Rinse and prime through the membrane to the filtrate outlet:

- PF 1000 N with approximately 500 ml
- PF 2000 N with approximately 1500 ml

Add 5000 IU heparin to the remaining 1000 ml of priming solution. Rinse with another 500 ml of the heparinized priming solution through the filtrate outlet.

Clamp the filtrate line and unclamp the venous line again. Rinse and prime the blood circuit with the remaining heparinized priming solution. Save 250 ml of the priming solution to rinse the blood circuit just before initiation of the treatment.

The plasmafilter is now ready for use.

TREATMENT

CAUTION !

In order to avoid hemolysis during plasmfiltration treatment:

- The pressure gradient between arterial inlet and filtrate outlet should be strictly controlled. **Do not exceed** a transmembrane pressure of:
 - PF 1000 N = 200 mmHg
 - PF 2000 N = 120 mmHg

PF 1000N:

The bloodflow rate should not exceed 200 ml/min or fall below 50 ml/min. The optimal range is between 80 - 130 ml/min.

PF 2000N:

The bloodflow rate should not exceed 250 ml/min or fall below 100 ml/min. The optimal range is between 150 - 200 ml/min.

- The venous hematocrit shall always be less than 65%.
- The plasmafilter should be in the vertical position with the arterial end at the bottom.
- If for any reason the initiation of the treatment is delayed, once the final rinsing of the plasmafilter has been completed, it is recommended that at least 250 ml of fresh sterile physiological solution is let through the plasmafilter.

If applicable, it is recommended to administer a loading dose of heparin to the patient 2 to 5 minutes before initiation of the treatment session.

Connect the arterial bloodline to the vascular access device.

Unclamp the lines and set the speed of the blood pump in accordance with the procedures utilized at the clinic. Discard as much as possible of the physiological solution from the system.

WARNING ! Do not leave the patient unattended during flush out of the rinsing fluid.

Stop the blood pump, clamp the lines and connect the venous line to the access of the patient.

Unclamp the lines, start the blood pump and set treatment variables in accordance with the needs and tolerance of the patient.

CAUTION !

- Continuous monitoring of all relevant parameters shall be maintained throughout the treatment session.
- To minimize the risk of clotting in the plasmafilter, never stop the blood pump before the filtration pump. If for some reason it is necessary to recirculate the extracorporeal blood through the plasmafilter, both the infusion and filtration pump must be turned off.

DISCONTINUING TREATMENT

RETURNING THE BLOOD TO THE PATIENT

Stop the infusion and filtrate pump.

Stop the blood pump.

Clamp the arterial blood line and the arterial needle. Remove the arterial line from the patient.

Connect the arterial line to a container with physiological solution.

Unclamp the arterial line. Start the blood pump and set the desired rinseback flow rate to give as much blood back as possible to the patient.

WARNING !

- **Do not** leave the patient unattended while the blood is being returned.
- **Do not** deactivate the air detection system on the machine until the blood in the extracorporeal circuit has been completely returned and the blood pump has been turned off. Air entering the patient may cause fatal air embolism.

Clamp the venous line repeatedly, to facilitate adequate emptying of the plasmafilter.

When the required amount of blood has been returned to the patient, stop the blood pump, clamp the venous line and the fistula needle. Remove the venous line from the patient.

General information

INDICATIONS

The plasmafilter is intended for use in therapeutic plasmfiltration, thus in diseases where removal of plasma components is indicated.

The plasmafilter can also be used in plasmfiltration procedures where filtered plasma is reinfused.

The size, weight, cardiac status and general physical condition of the patient must be evaluated by the prescribing physician before each treatment. The choice of the appropriate plasmafilter and associated equipment as well as the treatment operating parameters are the sole responsibility of the physician. Special attention must be paid in connection with pediatric use.

It is recommended that the plasmafilter is used together with a machine that complies with national or international standards.

CONTRAINDICATIONS

There are no absolute contraindications for use of the plasmafilter. The device should be used only on the direction of a physician who has evaluated all of the pertinent features of this device in relation to the individual patient.

Patients with bleeding tendencies must be closely supervised during the treatment. Treatment should be discontinued if acute hemorrhage that cannot be corrected (hematemesis, hemoptysis and melena) occurs during plasmfiltration.

ANTICOAGULATION

Anticoagulation may be necessary to prevent thrombus formation and clotting in the extracorporeal circuit. Anticoagulation therapy should be administered and monitored under the directions of a physician.

PERFORMANCE

The performance values given below should be regarded as approximate. When the plasmafilters are used under clinical conditions, different values may be obtained due to the clinical setting and measuring technique as well as to variations in the Polypropylene membranes and the manufacturing.

Type	PF 1000 N	PF 2000 N
Filtration in vitro (ml/min) (1) at TMP = 80 mmHg Q_b (ml/min)		
50	15 - 19	-
100	29 - 36	39 - 48
150	43 - 53	57 - 69
200	54 - 67	74 - 91
Sieving coefficient (total protein)	min. 0.95	min. 0.95
Flow resistance in vitro (mmHg) (1) Blood compartment at $Q_b = 100$ ml/min and TMP = 80 mmHg	max. 30	max. 35
Priming volume in vitro (ml) Blood compartment	23	41
Residual blood volume (ml)	< 1	< 1
Maximum TMP (mmHg) (TMP = $P_a - P_f$)*	200	120
Minimum blood flow (ml/min)	50	100
Maximum blood flow (ml/min)	200	250

acc. to DIN 58 353

(1): measured with bovine blood, hematocrit 32 %, protein content 60 g/l at 37° C. Upon request, information about methods used to obtain performance data is available.

* see Appendix: Pressure monitoring

SPECIFICATIONS

The blood and filtrate ports are designed and intended for use with Gambro blood lines. acc. to ISO 8637

Type	PF 1000 N	PF 2000 N
Effective membrane area (m²)	0.15	0.35
Fiber dimensions		
Inner diameter (µm)	330	330
Wall thickness (µm)	150	150
Max. pore size (µm)	0.5	0.5
Sterilization agent	Ethylene oxide	
Components	Materials	
Membrane	Polypropylene	(PP)
Potting material	Polyurethane	(PUR)
Housing and caps	Polycarbonate	(PC)
Protective caps blood side	Polypropylene	(PP)

* APPENDIX PRESSURE MONITORING

Normal TMP calculations are not relevant in plasmfiltration because there is a large pressure drop in the filter. This means that the pressure in the arterial end (bottom) of the filter is higher than the mean filter pressure. Since the TMP is normally calculated from the mean filter pressure this is lower than the arterial TMP.

P_a = Arterial pressure

P_v = Venous pressure

P_f = Filtrate pressure

The pressure gradient across the membrane at the bottom of the filter is defined as:

Arterial TMP = $P_a - P_f$

This should not be confused with the average transmembrane pressure which is defined as:

$$TMP = \frac{P_a + P_v}{2} - P_f$$

Warning: In order to minimize the risk of hemolysis do not exceed an arterial TMP for
PF 1000 N = 200 mmHg
PF 2000 N = 120 mmHg.

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