

Flame retardants

Application 2425 GC

Analysis of brominated diphenylethers and brominated biphenyls (flame retardants) via GC/MS following ROHS guideline

Technique : GC-Capillary
 Column : Fused Silica factorFour15m x 0.25mm, VF-5ms, df = 0.10 µm
 Part no: CP8938
 GC : CP3800
 Temperature : 100 °C, 2 min → 250 °C, 5 °C/min, 2 min → 300 °C, 5 °C/min, 2 min. → 320 °C, 5 °C/min, 6 min
 Carrier gas : He 0.8 mL/min, constant flow
 Injection : PTV 1078; Splitless, Temperature: 150 °C, 0 min → 350 °C, 200 °C/min
 Sample size : 2 µl, 20-2000 pg on the column
 Detector : MS-Saturn in ms/ms and MRM, Trap: 240 °C; Transferline: 285 °C;
 Solvent : isooctane

Courtesy : Rosemarie Stäbler, Dekra Umwelt GmbH, Stuttgart, Germany

Peak Identification

1	PBB2	2-bromobiphenyl	11	PCB209	internal Standard
2	BDE2	3-bromodiphenylether	12	PBB153	2,2',4,4',5,5'-hexabromobiphenyl
3	BDE15	4,4'-dibromodiphenylether	13	BDE153	2,2',4,4',5,5'-hexabromodiphenylether
4	PBB15	4,4'-dibromobiphenyl	14	BDE183	2,2',3,4,4',5,6-heptabromodiphenylether
5	PBB30	2,4,6-tribromobiphenyl	15	PBB209	decabromobiphenyl
6	BDE30	2,4,6-tribromodiphenylether	16	BDE203	2,2',3,4,4',5,5',6-octabromodiphenylether
7	PBB52	2,2',5,5'-tetrabromobiphenyl	17	PBB250	octabromobiphenyl (isomer)
8	BDE62	2,3,4,6-tetrabromodiphenylether	18	BDE206	2,2',3,3',4,4',5,5' 6-nonabromodiphenylether
9	PBB103	2,2',4,5',6-pentabromobiphenyl	19	BDE209	decabromodiphenylether
10	BDE99	2,2',4,4',5-pentabromodiphenylether			

