

# **Cerity Networked Data System for Pharmaceutical QA/QC**

**Revision A.02.02**

## **Specifications**

**July 2004**

### **Introduction**

Cerity NDS for Pharmaceutical QA/QC is a high availability networked data system specifically designed for pharmaceutical QA/QC laboratories. It is powerful, flexible and utilizes an integrated client-server architecture enabling seamless industry standard distributed client-server scalability. Its user interface is optimized to model the way analysts work in the QA/QC environment, fully supporting their everyday tasks. Cerity NDS for Pharmaceutical QA/QC is available in two configurations – A stand-alone system as entry level and a distributed client/server system.



**Agilent Technologies**

## Table of Contents

---

<b>Overview</b>	3
Certy NDS for Pharmaceutical QA/QC Professional	3
Certy NDS for Pharmaceutical QA/QC Client/Server	3
<b>General Description –</b>	
<b>Certy NDS for Pharmaceutical QA/QC Client/Server</b>	4
<b>Data Acquisition and Instrument Control</b>	4
Acquisition Controllers and Reprocessing Servers	4
Agilent Instrument Driver	4
Drivers for Non-Agilent Instruments	4
<b>Data Processing and Review</b>	5
<b>Software Administration Console</b>	5
<b>Network Infrastructure</b>	5
Overview	5
Recommended Topology	5
Recommended Network Settings	5
Network Interfaces	6
Assignment of IP Addresses	6
<b>Support for High Availability Configurations (Failure Resilience)</b>	7
Support for Clustering (Server Failover)	7
Certy Database Server	7
Acquisition Controller	7
<b>Supported Analytical Instrumentation</b>	8
Agilent 1100 Series Liquid Chromatograph	8
Agilent 1090 Series	8
Agilent 1090 Solvent Delivery Systems	8
Agilent 1090 Sampling Systems	8
Detector systems for the Agilent 1090	8
Interfacing of the Agilent 1090	8
Number of Instruments	8
Agilent 35900E Dual Channel Interface	8
Agilent 6850 and 6890 Gas Chromatographs	9
HP/Agilent 5890 GC	9
Supported Sampler Configurations	10
Water Alliance	10
Shimadzu LC 10A Series	10
<b>Certy Generic Instrument Control (GIC)</b>	11
<b>Hardware Requirements</b>	12
Certy NDS for Pharmaceutical QA/QC Professional	12
Certy NDS for Pharmaceutical QA/QC Client/Server	12
Database Server Configurations Overview	12
Terminal Server Configurations	12
Acquisition Controller	13
Review Client	14
<b>Operating System/Software Requirements</b>	15
Oracle Licensing	15
Certy NDS for Pharmaceutical QA/QC License	15
GMP Module License	16
Instrument Control and Spectral Evaluation Licenses	16
Installation Qualification Tool (IQT)	16
Operation Qualification Tool (OQT)	16
<b>Functional Specifications – Application</b>	18
<b>Functional Specifications – Electronic Records and</b>	
<b>Electronic Signatures Checklist (21 CFR Part 11)</b>	33

## Overview


### Cerity NDS for Pharmaceutical QA/QC Professional

Cerity NDS Professional is the solution for laboratories that require instrument control, data acquisition, data analysis, flexible reporting and support for up to eight 2D chromatographs (up to four 3D chromatographs) controlled by a single computer with strict adherence 21 CFR Part 11 (electronic records and electronic signatures) and related predicate rules such as to 21 CFR 210 (GMP) and 21 CFR Part 211 (cGMP). This configuration provides system access for one user at a time. It is designed for small laboratories that require secure data storage, data collection of electronic analytical records and support for several instruments without the need for multi-user operation and centralized chromatography data management.


The underlying technical infrastructure of Cerity Professional (e.g. processes, services and user interface components) is identical to the client/server configuration. A Cerity Professional system can be reused as a client when scaling up to a Cerity client/server configuration.

### Cerity NDS for Pharmaceutical QA/QC Client/Server

The Client/Server configuration extends the capabilities of Cerity Professional by distributing system tasks between the central database server, any number of connected instruments, acquisition controllers and review client workstations. This allows for multiple users to access and concurrently work with the central database and any connected



**Figure 1**  
Cerity NDS for Pharmaceutical QA/QC Professional allows a single user to control, acquire and process data from up to 8 dual channel instruments.



**Figure 2**  
Example configuration of Cerity NDS for Pharmaceutical QA/QC Client/Server with one central Oracle database server and two groups of instruments controlled by dedicated acquisition controllers. Users access the system through the review client PCs.

instrument transparently from any connected client workstation. You can configure as many clients

as necessary on the Cerity NDS system, considering any licensing and resource limitations.

## General Description – Cerity NDS for Pharmaceutical QA/QC Client/Server

---

In a client/server environment, the Cerity NDS database server is a Microsoft Windows server hosting an Oracle database. This is the central data repository of the Cerity NDS system. The system includes the following networked components: database server, acquisition controllers, instruments, printers, client workstations and other devices. All raw data (acquired sig-

nals), meta-data (i.e. methods, calibration information, calculation formulae, logbook), and calculated results are stored centrally in the database along with the computer-generated audit trail and revision information.

Standard queries allow searching, retrieving and displaying data for review and other purposes, such as

instrument status checking, results review, results approval and reporting.

The size of the database depends on the number of concurrent users, concurrent instruments and the amount of data available online in the database.

## Data Acquisition and Instrument Control

---

### Acquisition Controllers and Reprocessing Servers

Client/Server configurations allow the addition of optional acquisition controllers to balance instrument data stream and buffer it before upload to the central database server. The Cerity NDS acquisition controller is a dedicated Microsoft Windows Workstation running the Cerity NDS acquisition controller software. The acquisition controller performs data acquisition and instrument control. This component controls instruments that have been scheduled at the review client to execute an analysis using the specified method. It collects and processes the raw data and transmits it to the central database server. Acquisition controllers can theoretically be used as review clients. However, the background processing load on the acquisition controller may decrease system performance for interactive use and is therefore not recommended for routine use.

### Agilent Instrument Drivers

The suite of drivers to control the supported instruments is installed

with the Cerity NDS base software. The supported instruments are the Agilent 1100 HPLC, the Agilent 5890, 6850 and 6890 Gas Chromatographs, the Agilent 35900E A/D, the Shimadzu LC10vp LC and the Waters Alliance 2690/2695 LC with a Waters 2487 Dual Wavelength Detector.

Agilent instruments implement Level 4 instrument control using standard LAN communications (TCP/IP). Cerity NDS controls instrument parameters and components and collects digital signals from detectors. Data can be acquired at rates up to 200Hz, the data rate required for fast GC and supported on Agilent GCs.

One Agilent instrument control license is required for each instrument.

### Drivers for Non-Agilent Instruments

Cerity NDS provides full instrument control of and data acquisition from the Waters Alliance 2690 and 2695 LC system. Instrument control of and data acquisition

from the Waters Alliance LC require a direct connection between the computer and instrument, using an Agilent 82350 PCI high-performance GPIB interface adapter and cable. The standard interface control library software (SICL) required for GPIB communication is delivered with the Cerity NDS software. It must be installed on the acquisition controllers where the GPIB support is required.

A separate instrument control license (Agilent product number G4062AA) is required for each Waters Alliance HPLC system controlled by Cerity NDS for Pharmaceutical QA/QC.

Cerity NDS provides full instrument control of and data acquisition from the Shimadzu LC10vp HPLC, through Cerity's Generic Instrument Control interface in conjunction with an Agilent 35900E dual channel interface. Instrument control is performed through the serial port, and data acquisition through the analog/digital IO interface of the 35900E.

## Data Processing and Review

---

The Cerity NDS review client is the interface to the Cerity NDS system. Users perform their work with the graphical user interface

(GUI) of the review client. Sample and sequence scheduling, instrument status monitoring, method setup, results review and

approval and other tasks can be performed from the Cerity NDS review client.

## Software Administration Console

---

The Cerity Software Administration module provides a so-called “snap-in” for the Microsoft Management Console (MMC). MMC is a console designed to integrate management tools and functions and to present a common visual

environment for management applications. Only operating system administrators can log on to the console. The Cerity Software Administration module permits administrators to set up system components, administer and track

software licenses and license consumption, application modules, instruments, logon permissions, user capabilities and roles, auditing and system wide settings for electronic sign-off and numeric formatting.

## Network Infrastructure

---

### Overview

The Cerity NDS infrastructure is based on operating system and common IT standards. For instance, Cerity NDS is standardized on communication through the TCP/IP protocol for all Cerity network nodes, including analytical instruments, computers, and servers.

### Recommended Topology

The Cerity system relies on the computer network infrastructure. In order to provide optimum performance and uptime, the network must meet design criteria in terms of available bandwidth, IP address assignment, name resolution and isolation of the lab subnet from the corporate network.

Cerity acquisition controllers capture the real-time data of networked instruments and transfer it to the central database server, which is typically installed in a data center.

Cerity acquisition controllers are designed to automatically recon-

nect to the required server resources (i.e. Cerity system services and the database), but to protect instrument analysis data from LAN/WAN failures that affect the connection to the database server in the data center, Cerity acquisition controllers must be installed on the same subnet as the instruments. For most environments, this actually requires to place the acquisition controllers at the departmental level (i.e. the lab) and not in the data center.

### Recommended Network Settings

Available bandwidth settings and communication modes can vary, depending on the networking infrastructure, switches, and the network interface cards used in computers and instruments. Due to the lack of an established industry standard, auto-negotiation (auto-synchronization) of network speed settings has been reported to cause communication failures in some cases, especially where both computers and switches were configured to use this mode.

Agilent Technologies recommends a two-tiered approach: For dedicated laboratory systems, fixed network speed settings should be used at both ends. In this case, IT departments adjust the infrastructure settings as needed, e.g. when equipment is moved. If this level of control is not available, then fixed network speed settings are recommended for instruments and clients with the switches configured to auto-negotiate.

1. Typically, the supported infrastructure for client/server configurations of Cerity requires a switched environment. Switches have become an inexpensive way to properly isolate the infrastructure segments from each other. Hubs should remain an exception in a properly designed network, but can be accommodated on a case-by-case basis. Servers are typically connected to the backbone through a main switch. The interface should be set to use the maximum speed (at least 100 MBit/s) full duplex.

2. It is strongly advised that acquisition controllers should reside on the same subnet as the instruments. Exceptions to this rule can be justified on a case-by-case basis and may require additional measures for the infrastructure between the laboratory and the server. The recommended interface setting is at least 100 MBit/s full duplex.
3. The recommended network speed setting for PC clients is at least 100 MBit/s full duplex fixed at both ends.
4. The network speed supported by analytical instrumentation from Agilent depends on the exact model. The recommended setting for the Agilent 1100 HPLC, 6850 GC and the 35900E Dual Channel Interface, is 100 MBit/s half duplex. The Agilent 6890N Gas Chromatograph, is only supported at 10 MBit/s half duplex.
5. A simple recommendation is to set the network speed for all analytical instruments to 10 MBit/s half duplex. This allows moving instruments in the lab.
6. Network ports on the switches need to match the speed settings described above.
7. If fixed speed settings are not in line with existing IT procedures and practices, switches may use "autonegotiate".
8. Agilent Technologies strongly discourages from using auto-negotiate at both ends due to implementation incompatibilities observed with equipment from some manufacturers.

Agilent P/N LC	Agilent P/N GC	Description	Minimum FW
N/A	N/A	HP J2552B MIO card for 10Base-T, 10Base2	A.08.32
G1846A	G1847A	HP JetDirect 400N (MIO) card for 10Base-T, 10Base2, 100Base-TX	K.08.32
G1369A	N/A	Agilent LAN Communication card for the Agilent 1100 and 8453 Series. Adherence of G1369A to the requirements of IEEE 802.3, 2002 (Ethernet) has been verified by IOL, the Interoperability Laboratory of the University of New Hampshire's Research Computing Center.	All revisions.

**Table 1**  
**Supported instrument network interface cards.**

### Network Interfaces

To connect Agilent Technologies instruments to the Cerity NDS for Pharmaceutical QA/QC system, a suitable network interface card is required. Table 1 lists the network interface cards supported with the Cerity NDS for Pharmaceutical QA/QC software.

### Assignment of IP Addresses

- The system requires a static IP address for the database server module, which is installed by selecting the Cerity Professional or Cerity Database Server installation.
- Proper name resolution must be enforced for all Cerity devices to guarantee connectivity between all networked components of the system even if IP addresses get reassigned. This allows referencing network nodes by hostname rather than IP address.
- Agilent recommends that fixed IP addresses for clients, acquisition controllers and instruments be allocated.
- If fixed IP address assignment is not possible, Cerity clients and instruments may receive their IP addresses by using the BOOTP, or dynamically using the Dynamic Host Configuration Protocol DHCP.

- The Cerity NDS system includes a dedicated BOOTP server (Bootstrap Protocol) that operates as a service in Windows. The Agilent BOOTP service can be used if no other BOOTP server or DHCP server are available to assign IP addresses to networked instrumentation and clients.
- The instrument configuration requires specific buffer settings on the Jetdirect interface card. Software. The Cerity Software automatically configures these buffer settings.
- We recommend reserving IP addresses for Cerity devices and instruments when using DHCP. This ensures that a device will always have the same IP address. IP address reservation means that a static IP address is assigned to a dedicated MAC address.
- The Agilent 1100 HPLC and the Agilent 6890/6850 GC allow setting instrument IP addresses from the handheld controller or local keyboard. BootP or DHCP servers are not required if fixed IP addresses are assigned to the instrument via the local keyboard or user interface.



## Support for High Availability Configurations (Failure Resilience)

---

### Support for Clustering (Server Failover)

Cerity for Pharmaceutical QA/QC supports clustering in the so-called “Active-Passive Mode”, which is also referred to as “Active-Standby”. All Cerity permanent data is held on the shared storage in addition to the so-called quorum information needed by Microsoft clustering.

Microsoft Windows 2000 Advanced Server provides clustering at the operating system level. Microsoft Cluster Server (MSCS) runs services required for cluster operation. In an active-passive cluster, two servers (nodes) share access to an external storage. Both computers can see the data in the shared storage, but only one of the nodes is active while the other is effectively on standby, waiting to take over if the primary (and active) computer fails to respond. Both servers have physical IP addresses as well as hostnames by which they can be reached through the computer network. In addition, a separate virtual IP address and machine name are configured for the cluster. The virtual address is used to access the active physical node of the cluster. Required system services run on the active node during normal operation. In the case of a failover, those services are started on the inactive node.

Oracle Fail Safe (OFS) services make the Oracle RDBMS cluster-aware and interact with the operating system to respond to failover system events. OFS

ensures data integrity on the permanent data store in a failover scenario. OFS is an optional utility available for Oracle 9i. Cerity A.02.01 SR1 and higher support both MSCS and OFS. Oracle Enterprise Edition-specific functionality such as the Transparent Application Failover (TAF) and optional packages such as Real Application Clustering (RAC) provide additional support to handle in-flight data or database replication. These features are not supported by the Cerity NDS software.

### Cerity Database Server

Cerity system services do not interact directly with the operating system to manage clustering related events. Instead, Cerity respond to cluster failover scenarios through MSCS. Cerity services running on acquisition controllers and clients automatically reconnect to the database server cluster after the failover is completed. Typical failover delays have been measured at less than one minute. Cerity acquisition buffering functions may be triggered during the failover delay.

A working failover cluster configuration requires that all Oracle data files as well as the Cerity report storage be located on the shared storage of the cluster server. Increased hardware resources should be specified as per Microsoft recommendations. Specifically, a separate mirrored disk (a so-called quorum disk) is highly recommended.

### Acquisition Controller

Cerity supports enhanced acquisition buffering for instrument measurement data. In the event of database connectivity loss, acquired data is buffered on the acquisition controllers, and automatically spooled to the database once the database connectivity is re-established. To minimize buffering requirements, data analysis is suspended in acquisition buffering mode. Data analysis results can be generated by reprocessing the buffered raw data. The Cerity system continues and completes all sequences running at the time of the network failure.

## Supported Analytical Instrumentation

### Agilent 1100 Series Liquid Chromatograph

Full (Level 4) instrument control of the Agilent 1100 liquid chromatograph via LAN interface. In revision A.02.02 and higher, the following modules are supported:

- Agilent 1100 isocratic, binary and quaternary pumps
- Agilent 1100 VWD, MWD and DAD (in 2D-chromatography and 3D spectral acquisition mode mode)
- Agilent 1100 FLD (fluorescence detector)
- Agilent 1100 RID (refractive index detector)
- Agilent 1100 standard and thermostatted autosampler
- Agilent 1100 thermostatted column compartment
- Agilent 1100 vacuum degasser
- Agilent 1100 well-plate sampler and thermostatted well-plate sampler

One Agilent 1100 LC instrument control license is required per Agilent 1100 LC system connected to the Cerity software.

*Note: A LAN interface is required for the Agilent 1100 system. Cerity does not support GPIB for this instrument.*

### Agilent 1090 Series

- Agilent 1090 Series II/L
- Agilent 1090 Series II/M

*Note: The 1090A series is not supported.*

### Agilent 1090 Solvent Delivery Systems

- DR-5 pump system binary and ternary
- PV-5 pump system ternary and quaternary

*Note: The isocratic system will not be supported.*

Agilent P/N	Agilent 1100 Module Description	Minimum Firmware
G1310A	Isocratic pump	A.05.06
G1311A	Quaternary pump	A.05.06
G1312A	Binary pump	A.05.06
G1313A	Autosampler	A.05.06
G1329A	Thermostatted autosampler	A.05.06
G1330B	Autosampler cooling module	N/A
G1367A	Agilent 1100 well-plate sampler	A.05.07
G1368A	Agilent 1100 thermostatted well-plate sampler	A.05.07
G1314A	Variable wavelength detector (VWD)	A.05.06
G1315A	Diode-array detector (DAD)	A.05.06
G1315B	Diode-array detector (DAD)	A.05.06
G1365A	Multiple wavelength detector (MWD)	A.05.06
G1365B	Multiple wavelength detector (MWD)	A.05.06
G1321A	Fluorescence detector (FLD)	A.05.06
G1362A/B	Refractive index detector (RID)	A.05.06
G1316A	Thermostatted column compartment (TCC)	A.05.09
G1322A	Online degasser	N/A
G1323B	Control module	B.03.01

**Table 2**

**Supported Agilent 1100 Series LC modules.**

### Agilent 1090 Sampling Systems

- Standard autosampler with vial positions 0 to 99 and special position 100
- 25- $\mu$ L and 250- $\mu$ L syringe options

*Note: The cooled autosampler and manual injector will not be supported.*

### Detector systems for the Agilent 1090

- 1090 diode-array detector in 2D mode with up to 5 signals
- Other detectors supported by the Cerity software (e.g. 1100 Series detector)
- 35900E

*Note: The programmable filter-photometric detector will not be supported by Cerity if the 1090 is equipped with an FPD instead of a built-in DAD. This requires the user to disconnect the FPD.*

### Miscellaneous

- Cerity controls the four external contacts

### Interfacing of the Agilent 1090

- The 1090 will be interfaced via GP-IB to the acquisition controller. We support the 82350A and 82350B cards with firmware revision L.02.01 or higher.
- Agilent does not support Waters Alliance and Agilent 1090 to be configured on the same acquisition controller.

### Number of Instruments

- Cerity supports up to four 1090 instruments in 2D mode or up to two instruments in 3D mode connected to a single acquisition controller.

### Agilent 35900E Dual Channel Interface

This interface can be used to acquire up to two independent channels of data from instrumentation that is not directly controlled by the software. It supports data rates up to 100 Hz. The Agilent 35900E Dual Channel Interface



supports BCD-coded vial number and can be used to track the vial position of each injection from third party auto samplers in the software. One Agilent 35900E A/D instrument control license is required per Agilent 35900E dual channel interface connected to the Cerity software.

A LAN interface is required for the dual channel interface. The GPIB interface is not supported in Cerity NDS for Pharmaceutical QA/QC for this device, refer to table 3 for firmware requirements.

#### Agilent 6850 and 6890 Gas Chromatographs

Cerity NDS for Pharmaceutical QA/QC supports the Agilent 6890A, and 6890N (table 4).

A LAN interface is required for the Agilent GC. The GPIB interface is not supported in Cerity NDS for Pharmaceutical QA/QC for this instrument (table 5).

Cerity NDS for Pharmaceutical QA/QC supports the 6850 GC hardware listed in table 6.

#### HP/Agilent 5890 GC

Revision A.02.02 of Cerity for Pharmaceutical QA/QC implements a framework for generic instrument control. This framework is used to provide instrument control and data acquisition with the Agilent 5890 GC.

- HP 5890 Series II and HP 5890 Plus
- The 19257–60020 card (revision C or higher) is needed for controlling the instrument
- GC ROM firmware revision “HP5890A.03.00” or higher

P/N	Description	Firmware	Comment
35900E	A/D converter	E.01.02	Data acquisition, remote start/stop, BCD

**Table 3**  
Minimum firmware requirements .

Agilent 6890A GC			
Inlet:	EPC—S/S, EPC—P/P, EPC—COC		
Column inlet:	Front, back, unspecified		
Column outlet:	Front, back, other, MSD, AED		
Oven:	High ramp rate		
Detector:	AIB, EPC—FID, EPC-FPD, EPC-NPD, EPC—TCD, EPC—ECD, EPC—uECD		
Cryo:	CO <sub>2</sub> , N <sub>2</sub>		
Valves:	GSV, LSV, Multiposition, Switching		
Aux:	EPC, Temperature		
Data channels:	In revision A.02.02, dual simultaneous (“dual tower”) injection is not supported		

**Table 4**  
Supported 6890 GC configurations.

Agilent P/N	Description	Firmware	Comment
Gas Chromatographs			
G1530A	6890A	A.03.05	With Jetdirect card
G1530N	6890N	N.04.05	Network ready 6890 GC
G2630A	6850	A.03.01	With Jetdirect card
G2630A	6850	A.05.01	Network ready 6850 GC
Automatic Liquid Samplers			
G1512A	7673 Controller	A.01.09	This controller is required to use 7683B injector and tray with the for 6890A GC. 6890N: firmware N.05.02 or greater required. 6850N: firmware A.05.02 or greater required. The 7683B autoinjector has two additional parameters: Wash mode and solvent saver mode. By default, these parameters are turned off. To use them, they have to be configured on the GC keyboard before connecting the GC to Cerity for Pharmaceutical QA/QC. 6890N: firmware N.05.02 or greater required. 6850N: firmware A.05.02 or greater required.
G1513A	7673 Injector	A.10.04	
G1514A	7673 Tray	N.A.	
G2612A	7683A Controller	A.01.07	
G2613A	7683A Injector	A.10.04	
G2614A	7683A Tray	A.01.01	
G2912A	7683B Controller	A.02.00	
G2913A	7683B Injector	A.11.00	
G2916A	7683B Tray	A.02.00	
18593B	7673 Injector	A.08.00	
18596C	7673 Tray	n/a	
Headspace Samplers			
G1289B	7694B	1.04	
G1290B	7694B	1.04	
G1888A	G1888A	A.01.02	
	Headspace sampler		

**Table 5**  
Agilent GC minimum firmware requirements.

- Inlets Split/Splitless and Cooled On Column. The inlets are supported with or without EPC.
- Electronic Pressure Control: 2-Channel or 6-Channel EPC  
The 5890 VSIA supports EPC Channels A, B, C, D, E and F.
- Detectors: TCD, FID, NPD, ECD and FPD. Time-programmable polarity switching for the TCD is not supported.
- Oven: Standard oven with or without cryogenic cooling

#### Supported Sampler Configurations

- ALS controller G1512A firmware revision A.01.05 or higher
- ALS injector 18593B and G1513A 8 Vial high-density turret is not supported
- Sample tray 100-vial tray (18596B/C) barcode reader is not supported
- The new 7683B injector (G2913A) does not provide a serial interface and is not supported with this driver

#### Waters Alliance

Full control of the solvent delivery system, column heater and autosampler of the Waters 2690 and 2695 Alliance Liquid Chromatograph via GPIB interface (optionally available with the software). Up to four Waters Alliance chromatographs can be controlled through a single GPIB interface. Full control of the Waters 2487 Dual Wavelength Detector via GPIB interface. One instrument control license is required per Waters Alliance system (consisting of one Waters 2690/2695 Alliance mainframe and one Waters 2487 detector) connected to the Cerity software.

#### Agilent 6850 GC

Injector:	G2613A/G2880A
Inlet:	Split/splitless, purge/pack
Column Inlet:	Inlet, unspecified
Column Outlet:	Detector, MSD
Detector:	FID, TCD
Cryo:	CO <sub>2</sub> , N <sub>2</sub>
Valves:	GSV, LSV, multiposition, switching
Aux:	Temperature
Data Channels:	1
Handheld controller	G2629A

**Table 6**  
**6850 GC hardware supported by Cerity NDS for Pharmaceutical QA/QC.**

Model Number	Description	Firmware	Comment
Waters 2690	Waters Alliance	2.02 or higher	2.02 is required to run Cerity compliance for the Waters Alliance. 2.02 is required to run Cerity compliance for the Waters Alliance.
Waters 2695	Waters Alliance	2.0 or higher	
Waters 2487	Waters Dual Wavelength Detector	1.03	

**Table 7a**  
**Waters instrument minimum firmware requirements.**

Module	Description	Firmware
SIL-10ADvp	Autosampler	5.32
LC-10ADvp	Pump	5.26
FCV-10ALvp	Switching valve for quat. gradient	n.a.
DGU-14Avp	Degasser	n.a.
SCL-10Avp	System controller	5.32
SPD-10AVvp	UV/VIS detector	5.24
SPD-10Avp	UV/VIS detector (identical to SPD-10AVvp but without a VIS lamp)	5.24
CTO-10ACvp	Column compartment	5.24
CTO-10ASvp	Column compartment	5.24

**Table 7b**  
**Shimadzu LC 10A Series firmware requirements**

A remote start-stop cable is required to synchronize the Waters detector with the LC. See table 7a for firmware requirements.

#### Shimadzu LC 10A Series

Cerity controls the Shimadzu LC 10A series through the Cerity Generic Instrument Control framework (GIC).

## Cerity Generic Instrument Control (GIC)

---

Agilent's Generic Instrument Control (GIC) concept is a new approach to standardize generic instrument control for non-Agilent instruments. GIC allows to integrate non-Agilent chromatographic instruments into the Cerity architecture. GIC is integrated into the core of the Cerity software and is always installed along with Cerity. Through the GIC framework, Cerity sees and controls a third party instrument similar to a specific known instrument type such as the 1100 Series LC. The GIC concept allows to develop drivers for third party instruments without changing the Cerity source code. Device drivers, called Vendor Specific Instrument Adapters (VSIA), can be developed independently of Cerity, because the interaction between the device driver and the rest of the system is limited to a well-defined interface. This is a clear advantage for companies such as independent hardware vendors, who want to develop their own device drivers.

Cerity GIC consists of two main components: the GIC framework and the VSIA. The GIC framework is installed as part of Cerity for the Pharma Client/Server system software. GIC is responsible for handling all the dynamic and configurable user interfaces in the instrument and method context, which are dealing with instrument control, status and error handling and instrument setpoints as part of the Cerity method.

VSIA is the real “instrument driver”.

It is responsible for the communication to the instrument (downloading method setpoints, reading actual values, error handling etc.) and the translation of the set point values into the control commands the instrument understands. VSIA supports “plug-and-play” installation. “Plug-and-play” allows the installation of VSIA files on a running Cerity for Pharma Client/Server system without interrupting the operation. Every third party instrument needs its own specific VSIA. Every VSIA must be installed separately.

The communication between the GIC framework and the VSIA is done using the Common Instrument Control Language (CICL). CICL specifications combine a state model, a set of commands and a COM interface. Device drivers that adhere to the CICL specifications can be plugged into the Cerity software.

The GIC framework provides a ready-made library for VSIA to interpret GIC commands. All CICL commands are in XML format. XML provides great flexibility for communicating device specific parameters, and eases the writing of VSIA's responses. VSIA understands the CICL commands and translates them to a instrument specific protocol and vice versa.

The definition of all setpoints is done using an industry standard XML file format, which makes it easy to describe the instrument capabilities and setpoints using off-the-shelf tools. The instrument spe-

cific parameters are specified in a device configuration XML file called as VSIA XML. This is the starting point, where the GIC recognizes a new instrument type and its properties. The XML configuration file serves as a reference for the installation of the new instrument driver.

Sample sources are provided in the developer kit from Agilent to bring additional clarity. The communication between the VSIA and the instrument is established using the Agilent analog/digital converter AD 35900E interface as communication interface. In addition, the Agilent AD 35900E is responsible for acquiring the analog signal and converting it into digital data. The standard 10/100 LAN interface connects the Agilent AD 35900E to the data system. The implementation of VSIA into easy to understand generic standards allows advanced customers (future) or an experienced third party companies (today) to program the instrument control portion outside of Agilent.

The GIC framework covers the following device categories and makes them consistent with other Cerity drivers:

- pump
- injector
- column
- detector
- integrated instrument
- auxiliary devices

## Hardware Requirements

The data system consists of a personal computer (PC) and Agilent software. All hardware and peripherals must appear in the appropriate Microsoft Windows Compatibility Lists for the operating system. Cerity NDS for Pharmaceutical QA/QC is designed to run on computers that conform to the specifications listed below.

### Cerity NDS for Pharmaceutical QA/QC Professional

Table 8 lists the minimum computer configuration that is supported for a Cerity NDS for Pharmaceutical QA/QC professional system. This table provides a guideline for computer hardware specifications, such as the amount of random access memory (RAM), disk space measured in gigabytes (GB), central processing unit (CPU) speed measured in megahertz (MHz), etc.

### Cerity NDS for Pharmaceutical QA/QC Client/Server

#### Database Server Configurations Overview

Table 9 specifies the minimum and recommended hardware configurations for a database server.

#### NOTE:

- It is NOT recommended to use the Cerity database server as a print server for the Review Client computers.
- It is not recommended to use the Cerity database server for network administration services such as control.

#### Terminal Server Configurations

Table 10 specifies the minimum hardware requirements for distributing the Cerity NDS user

<b>Operating system type(s)</b>	<b>Microsoft Windows</b>	
<b>Operating system version(s)</b>	Windows 2000 Professional	Windows XP Professional
<b>Service pack version(s)</b>	Service Pack 4 MS Hotfix KB824146	Service Pack 1a MS HotFix KB824146 & KB823980
<b>Web browser</b>	Microsoft Internet Explorer 6.0 Service Pack 1	
<b>RDBMS</b>	Oracle 9i Client Release 2 Oracle 9i Patch Set 2 Oracle OLE DB 9.2.0.2 Patch	
<b>Configuration specifications</b>	<b>Low End</b> 4 single channel instruments, with a max. of 2 instruments that acquire on-line spectra, e.g. DAD or FLD	<b>High End</b> 8 single channel instruments, with a max. of 4 instruments that acquire on-line spectra, e.g. DAD or FLD
<b>Processor type or minimum speed, etc.</b>	Pentium 4, 1.2 GHz	Pentium 4, 2 GHz
<b>Minimum memory</b>	768 MB	1.5 GB
<b>Minimum mass storage</b>	80 GB	120GB
<b>Peripherals required</b>	Keyboard, mouse, CD-ROM drive	
<b>Interfaces</b>	Network interface	
<b>Display</b>	1024 x 768 pixels, 65536 colors	

**Table 8**  
Recommended hardware configuration for Cerity NDS Professional.

	Entry Level	Medium	High End
<b>Specification</b>	≤ 20 instrument channels ≤ 10 concurrent users ≤ 100.000 Tests	≤ 70 instrument channels ≤ 20 concurrent users ≤ 330.000 Tests	> 70 instrument channels > 20 concurrent users > 700.000 Tests
<b>CPU</b>	Xeon, 2.4 GHz	Xeon, 2.4 Dual processor recommended for >50 instruments	Xeon, 2x 2.4 GHz
<b>Memory</b>	1 GB	1.5 GB	3 GB
<b>Disk space</b>	3 GB Cerity server application and Oracle RDBMS installation 4,5 GB Data Table spaces for data, index and BLOBs 10 GB File share for archives and reports 50 GB Temporary space and archive redo log files	3 GB 140 GB 32 GB 150 GB	3 GB 290 GB 68 GB 300 GB
<b>Total (rounded)</b>	110 GB	330 GB	670 GB

**Table 9**  
Recommended hardware configurations for Cerity NDS database server.

<b>CPU</b>	<10 concurrent clients	<25 concurrent clients
<b>Computer type</b>	Dual Pentium 4, 1.4 GHz Desktop PC or server	Dual Xeon, 1.4 GHz Server
<b>Memory (RAM)</b>	768 MB	2 GB
<b>Disk size</b>	18.2 GB	18.2 GB ; second disk may be added for redundancy

**Table 10**  
Hardware prerequisites for Cerity NDS acquisition controller

interface through a dedicated Terminal Server. One or more Windows 2000 servers with Terminal Services and Citrix Metaframe XP can be configured to run the Cerity review client software. However, the terminal servers are not expected to run other software besides Cerity and cannot be configured in server farm. A Terminal Server configuration of Cerity for Pharmaceutical QA/QC will require approximately 200 MB of memory per user (RAM and/or virtual memory) and 6 3D chromatography instruments with spectral data acquisition.

### Acquisition Controller

In client/server configurations with multiple instruments, one or multiple separate acquisition controllers are required for instrument data acquisition and instrument control.

There is no hard-coded limit on the number of instrument channels that can be collected by a single acquisition controller. In fact, feasibility tests have shown that the system can control as many as 30 3D systems. However, the recommended maximum number of chromatography instruments for an acquisition controller is 15 2D chromatography instruments (i.e. without spectral data acquisition) and 6 3D chromatography instruments with spectral data acquisition.

Acquisition controllers can also be used to balance the reprocessing load within a distributed Cerity NDS installation. In this case, the acquisition controller serves as a dedicated reprocessing server.

	Low end	High end
<b>Processor type or minimum speed, etc.</b>	Pentium 4, 1.2 GHz	Pentium 4, 2.4 GHz
<b>Minimum memory</b>	1.5 GB	2.5GB
<b>Maximum number of instruments</b>	8 instruments with no on-line spectra capabilities or 3 instruments that acquire on-line spectra, like DAD or FLD	15 instruments with no on-line spectra capabilities or 6 instruments that acquire on-line spectra, like DAD or FLD
<b>Minimum mass storage</b>	40 GB	
<b>Peripherals required</b>	Keyboard, mouse, CD-ROM drive	
<b>Interfaces</b>	Network Interface	
<b>Display</b>	1024 x 768 pixels, 65536 colors	

**Table 11**  
**Hardware prerequisites for Cerity NDS review client.**

1. Plan your client/server systems to that you use no more than 15 data channels for each acquisition controller. For example, a 35900E with 2 data channels configured and a DAD with all five chromatography signals configured leaves eight available channels.
2. For large systems with a heavy processing load, add an additional acquisition controller for every 10 review clients in your client/server system to perform off-line reprocessing. Do not add any instruments to these systems. By default, the Cerity server assigns the reprocessing servers automatically. System administrators can explicitly assign a dedicated reprocessing server to specific review clients. This configuration is useful if certain review clients create a much higher data load than others.

The following table specifies the minimum hardware requirements for an acquisition controller.

If server computers are available for acquisition controllers, RAID controllers are recommended. Fast PCs with sufficient RAM and virtual memory are also acceptable. For system redundancy, please consider standby computers to be used in case of computer hardware failures.

To have Windows XP choose the best paging file size, choose "System managed size". The recommended minimum size is equivalent to 1.5 times the amount of RAM on your system, and 3 times that figure for the maximum size. Example, if you have 4 GB of RAM, the minimum paging file size is 6 GB and the maximum paging file size is 18 GB.

If the paging file reaches its maximum size in Windows 2000, a warning is displayed and the system may halt. To see whether your paging file is approaching its upper limit before it reaches the upper limit, check the actual file size and compare it to the maximum paging file size setting in the System utility in Control Panel. If these two numbers are close in value, consider increasing initial paging file size or running fewer programs.

### Review Client

Review clients are workstations used for interactive entry of sample and sequence data, method entry and management, scheduling of analyses, data review, reporting and result approval.

Table 12 specifies the minimum hardware requirements for review clients. Keep in mind, as previously stated, the acquisition controller and review client may operate together on a computer. It is recommended to deploy dedicated acquisition controllers for optimum performance and load balancing in the Cerity NDS system.

<b>Operating System type(s)</b>	Microsoft Windows	
<b>Operating System version(s)</b>	Windows 2000 Professional	Windows XP Professional
<b>Service pack version(s)</b>	Service Pack 4	Service Pack 1a
	MS Hotfix KB824146	MS HotFix KB824146 & KB823980
<b>Web Browser</b>	Microsoft Internet Explorer 6.0 Service Pack 1	
<b>RDMBS</b>	Oracle 9i Client Release 2	
	Oracle 9i Patch Set 2	
	Oracle OLE DB 9.2.0.2 Patch	
<b>Processor type or minimum speed, etc.</b>	Pentium III, 800 MHz	
<b>Minimum memory</b>	Minimum 256 MB, 512 MB recommended	
<b>Minimum mass storage</b>	20GB	
<b>Peripherals required</b>	Keyboard, mouse, CD-ROM drive	
<b>Interfaces</b>	Network interface	
<b>Display</b>	1024 x 768 pixels, 65536 colors	

**Table 12**  
**Hardware prerequisites for Cerity NDS review clients.**



## Operating System/Software Requirements

Table 13 defines the operating system requirements for each of the Cerity NDS for pharmaceutical QA/QC components. Table 14 defines the third party required software revisions to properly operate Cerity NDS for pharmaceutical QA/QC client/server and professional systems. Any of these components are installed automatically during Cerity setup, except for the Oracle Database Management System (DBMS) which needs to be installed separately.

### Oracle Licensing

Cerity NDS for Pharmaceutical QA/QC uses the Oracle RDBMS to manage and store its records.

- The Oracle RDBMS software may only be installed and used if the appropriate software licenses have been purchased. You must possess an Oracle license for each user account ("named user") established in your Agilent networked data system valid for use with the Agilent NDS software.
- The base products of Cerity NDS for Pharmaceutical QA/QC (Agilent G4000AA and G4001AA) include five application specific named user Oracle client licenses. These licenses are subject to a restricted use license and can only be used in conjunction with the NDS application.
- Agilent provides support for included Oracle software according to the application requirements of the respective Agilent networked data system. Further software maintenance for Oracle software must be purchased separately
- Alternatively, you may purchase full use Oracle licenses from

	Windows 2000 SP3 Professional	Windows 2000 SP3 Advanced Server	Microsoft Windows XP SP1/1a Professional
<b>Cerity NDS for Pharmaceutical QA/QC Database Server</b>	No	Yes	No
<b>Cerity NDS for Pharmaceutical QA/QC Review Client</b>	Yes	Yes	Yes
<b>Cerity NDS for Pharmaceutical QA/QC Acquisition Controllers</b>	Yes	Yes	No
<b>Cerity NDS for Pharmaceutical QA/QC Professional</b>	Yes	Yes	No

**Table 13**  
Cerity for Pharmaceutical QA/QC operating system requirements.

<sup>1</sup> Cerity NDS for Pharmaceutical QA/QC is fully tested and supported on the US-English, Japanese, French, German, Italian, Japanese and Spanish versions of Windows. On non-US English versions of the operating system, language-specific hotfixes of the operating system may be required. Any mandatory hotfixes are supplied on the Cerity NDS CD media. For correct interpretation of numeric floating point entries, specific options need to be configured in Windows' regional settings. Please consult the software status bulletin and release notes (readme.txt) for more information on regional settings, number and date formats required for correct handling of floating point number input.

Manufacturer	Product	Revision	Comment
Oracle	Oracle 9i DBMS	9.2.0.3.0	Shipped on product CD-ROM. Needs to be installed separately. Oracle Standard 9i is contained on Cerity product CD #2 through 6. Included with product.
Microsoft	Data Access Components (MDAC)	2.7	
Microsoft	Internet Explorer	6 SP1	Needs to be installed separately. Included on Cerity CD#1.
Agilent	SICL library	L.02.01	Required for the control of Waters Alliance Contained on product CD #1 .
Microsoft	Visual Basic	6.0 SP4	Runtime library
Microsoft	Visual C++	6.0 SP4	Runtime library

**Table 14**  
Other software requirements for Cerity for Pharmaceutical QA/QC.

- Oracle Co. or their authorized distribution partners.
- Each individual with a logon to the Cerity software requires a separate Oracle client license.
- Additional Oracle licenses, can be purchased using Agilent product number G1411A. Please note that Oracle licenses delivered by Agilent Technologies are application-specific and may only be used within the context of the Agilent networked data system.

### Cerity NDS for Pharmaceutical QA/QC License

The number of Cerity NDS concurrent licenses in use must not exceed the number of licenses installed, otherwise the license agreement is violated. Additional licenses are easy to order and easy to install. The Cerity NDS licenses float and are consumed by concurrent users. For the client/server configuration, there must be an equal number of Cerity NDS licenses as there are Cerity NDS concurrent users.

Cerity NDS for Pharmaceutical QA/QC Professional (Agilent G4000AA) is a single-user system and includes one concurrent Cerity user license. The following restrictions apply:

- Only one Cerity user can be logged on to the computer at any given time.
- G4002AA (add concurrent Cerity NDS user) is not applicable to the professional system. If the professional system is used by different operators at different times, each individual user requires a separate named user license of Oracle Standard valid for the Agilent NDS family of products.

#### **GMP Module License**

- One GMP license is needed for each concurrent user license (Agilent G4002AA)
- The GMP module enables strict auditing, audit comments and e-signature.
- The GMP module also enforces a strict results review/approval process. This will ensure that analysts review their own results before a peer reviewer and a final approval is given.
- This is an enabling license. No additional software is installed. The GMP module enables the audit node in the Cerity System Administration Console.

*Note: Only one GMP license is needed in the professional system.*

#### **Instrument Control and Spectral Evaluation Licenses**

One instrument control license is required per instrument controlled by the Cerity NDS software. The licenses are easy to install, and they are monitored by the application. Cerity instrument control

licenses are available for the following products:

- Cerity-P 1100 Series LC instrument control license (G4061AA)
- Cerity-P Waters Alliance LC instrument control license (G4062AA)
- Cerity-P 6890/6850 GC instrument control license (G4063AA)
- Cerity-P 35900E Dual Channel Interface instrument control license (G4064AA)
- Cerity-P Spectral Acquisition and Evaluation License (G4031AA). Enables data acquisition and data evaluation for spectral data obtained from an Agilent 1100 instrument with 3D spectral capabilities (diode-array detector or fluorescence detector). One spectral processing license is required per 3D instrument.
- Cerity-P Agilent 1090 HPLC instrument control license (G4067AA)
- Cerity-P Generic Instrument Control framework license (G4065AA) to enable instrument control and data acquisition through the Agilent 35900E dual channel interface using the Cerity Generic Instrument Control framework. Requires one license and one 35900E channel per GIC instrument. G4064AA license is included, but 35900E hardware is not. Requires one additional vendor specific instrument adapter for each instrument.
- Cerity-P Shimadzu LC10Avp VSIA instrument control license (G4067AA)
- Cerity-P Agilent 5890 VSIA instrument control license (G4068AA)

#### **Installation Qualification Tool (IQT)**

The Installation Qualification Tool is a computer-based qualification utility used to perform Installation Qualification (IQ) of the Cerity NDS for Pharmaceutical QA/QC system. Computer-based installation qualification protocols verify the completeness and intactness of the Agilent software installed on the PC. The computer-based installation qualification utility available for Cerity NDS for Pharmaceutical QA/QC reads required details from the system directly and inserts them into the document automatically. The utility provides input forms for details that cannot be extracted automatically from the system (software and hardware). The entry forms support further techniques for automated data entry such as bar coding. Execution of the software IQ protocol requires a valid IQ license. Without a valid license number, the final acceptance protocol cannot be generated.

#### **Operation Qualification Tool (OQT)**

OQ allows qualification tests at defined intervals on the data system and the connected instruments. Without a valid OQ/PV test result, the system must not be used. OQ/PV typically requires a series of different tests depending on the instrument, the lab's specification and the configured software capabilities. System IQ and OQ/PV are provided both as a product and as a service from Agilent Technologies. The scope of the validation services and products for Cerity NDS includes the Agilent 1100 HPLC, Agilent 6890/6850 GC, Agilent 35900E A/D and software system qualification. The Operational Qualification Tool is a computer-

based qualification utility used to perform Operational Qualification (OQ/PV) of the Cerity NDS for Pharmaceutical QA/QC system. The computer-based OQ protocols available for Cerity NDS for Pharmaceutical QA/QC uses well defined interfaces (so-called test harnesses) specifically designed into the software for the purpose of executing critical system test cases automatically. This comprises:

- Automatic low and mid-level functional tests that verify fundamental system-level functions that are not even covered by the traditional interactive protocols available for other data systems
- Automatic high-level system operation tests that verify application functionality such as sequencing, quantification or recalibration. These tests execute in unattended mode and the evaluation is performed automatically using known data source, prerecorded acceptance limits and self-evaluating reports.
- A number of test cases require scripted manual tests because of their interactive nature. The test scripts cover areas such as challenging logon security, auditing of interactive changes, authority checks, and archive/restore functions.

Execution of the OQ protocols requires a valid OQ license. Without a valid license number, the final acceptance protocol cannot be generated.

## Functional Specifications – Application

### Technology and Architecture

General description	Cerity NDS for Pharmaceutical QA/QC is a high availability, fully scaleable networked data system for analytical QA/QC laboratories that require chromatographic instrument control, data acquisition, data analysis, flexible reporting and with strict adherence 21 CFR Part 11 (electronic records and electronic signatures) and related predicate rules such as to 21 CFR 210 (GMP) and 21 CFR Part 211 (cGMP).														
Revision history	<p>Cerity for Pharmaceutical QA/QC has been available in the following revisions:</p> <table> <tr> <th>Revision number</th><th>Description</th></tr> <tr> <td>Revision A.01.01 (June 2002)</td><td>Initial release</td></tr> <tr> <td>Revision A.01.02 (August 2002)</td><td>Support for Windows 2000, support of 6850 GC</td></tr> <tr> <td>Revision A.01.03 (November 2002)</td><td>Support under Citrix Terminal server</td></tr> <tr> <td>Revision A.02.01 (March 2003)</td><td>Support for Oracle9i RDBMS and Windows XP for clients. Data model change to manage binary large objects entirely inside the Oracle database. Support for clustering of database server</td></tr> <tr> <td>Revision A.02.01 SR1 (June 2003)</td><td>Enable data buffering on acquisition controllers in case of network or server failure.</td></tr> <tr> <td>Revision A.02.02 (January 2004)</td><td>Enable 3D spectral data acquisition and data evaluation for diode array and fluorescence detectors</td></tr> </table>	Revision number	Description	Revision A.01.01 (June 2002)	Initial release	Revision A.01.02 (August 2002)	Support for Windows 2000, support of 6850 GC	Revision A.01.03 (November 2002)	Support under Citrix Terminal server	Revision A.02.01 (March 2003)	Support for Oracle9i RDBMS and Windows XP for clients. Data model change to manage binary large objects entirely inside the Oracle database. Support for clustering of database server	Revision A.02.01 SR1 (June 2003)	Enable data buffering on acquisition controllers in case of network or server failure.	Revision A.02.02 (January 2004)	Enable 3D spectral data acquisition and data evaluation for diode array and fluorescence detectors
Revision number	Description														
Revision A.01.01 (June 2002)	Initial release														
Revision A.01.02 (August 2002)	Support for Windows 2000, support of 6850 GC														
Revision A.01.03 (November 2002)	Support under Citrix Terminal server														
Revision A.02.01 (March 2003)	Support for Oracle9i RDBMS and Windows XP for clients. Data model change to manage binary large objects entirely inside the Oracle database. Support for clustering of database server														
Revision A.02.01 SR1 (June 2003)	Enable data buffering on acquisition controllers in case of network or server failure.														
Revision A.02.02 (January 2004)	Enable 3D spectral data acquisition and data evaluation for diode array and fluorescence detectors														
Description of data repository	Central, secure data repository based on Oracle database management system.														

### Supported Database Management System Oracle® 9i (9.2.0.3.0)

Data model	Object-relational data model.
Management of binary large objects (BLOBs)	Binary data (e.g. raw instrument data) is stored in so-called BLOBs and managed within the Oracle 9i database.
Characteristics of the distributed Cerity configurations	<ul style="list-style-type: none"> <li>• Full scalability (1 to n)</li> <li>• High availability configurations (acquisition buffering, server clustering)</li> <li>• Centralized system administration</li> <li>• Centralized backups</li> <li>• Import and export of data is possible between database servers</li> <li>• Multiple Cerity database servers may exist in a domain</li> <li>• Cerity review clients can connect to more than one Cerity servers</li> <li>• The number of review clients is limited by physical resources and available licenses only</li> <li>• Multiple acquisition controller are supported per Cerity database server</li> <li>• Up to 15 2D instruments or up to 6 3D instruments are supported per acquisition controller</li> </ul>
Maximum size of data base supported	No hard-coded limits. Three standard database configurations are offered (small, medium, large).
Design language/tools	UML, Visual C++, Visual Basic

### Operating System

- Microsoft Windows 2000 Advanced Server (Cerity C/S Database server)
  - Windows 2000 Professional Service Pack 4 (Professional, acquisition controller, review clients)
  - Windows XP Professional Service Pack 1a (Professional, acquisition controller, review clients)
- All configuration require Microsoft Internet Explorer 6.0 Service Pack 1
- For more details, refer to the respective chapter in this document.

continued ...

Compatibility with Windows Terminal Server	Client workstations of Cerity for Pharmaceutical QA/QC rev. A.01.03 and higher are supported for operation in a Windows 2000 Terminal Server environment. Thin client configurations are dictated by the Terminal Server provider (Citrix Metaframe or Microsoft) Agilent does not add any requirements unless specified otherwise. For details on the requirements and configurations of Cerity for Pharmaceutical QA/QC, please refer to a separate Technical Note, available as Agilent publication G4000-90100.
Other client/server capabilities	<ul style="list-style-type: none"> <li>• During running analyses, users can log out and log in without interruption of the sequence.</li> <li>• After log-out of the session, the user's license is released and available for another concurrent user.</li> </ul>
Hardware requirements	See separate chapter in this document.

### High Availability

	<p>The Cerity system supports the following failure resilience mechanisms:</p> <ul style="list-style-type: none"> <li>• Server failover using Microsoft Cluster Server and Oracle Failsafe configurations to ensure continued uptime even if the server fails</li> <li>• Acquisition buffering on the acquisition controllers ensuring that no analysis data can be lost even if server/ network infrastructure fails.</li> </ul>
Microsoft Windows 2000 Advanced Server Clustering	Supported in Active Standby mode
Oracle Failsafe (OFS)	Supported
Oracle Transparent Application Failover (TAF)	Not supported
Oracle Real Application Clustering (RAC)	Not supported

### Scalability

Cerity NDS for Pharmaceutical QA/QC Professional	One user at any one time.
Maximum number of concurrent users	Note: By definition, Cerity NDS for Pharmaceutical QA/QC Professional is a single-user, multi-instrument configuration.
Cerity NDS for Pharmaceutical QA/QC client/server	No hard-coded limits (function of database server configuration and available network bandwidth).
Maximum number of concurrent users	Typical configurations have 30-50 users and 80-100 instruments.

### User Interface

Characteristics of the user interface	<p>Graphical user interface designed in adherence to Windows standards and configurable to laboratory specific workflow (based on user roles and analysis specific requirements).</p> <p>The user interface of Cerity NDS for Pharmaceutical QA/QC is streamlined to adhere to the requirements of GMP regulated QA/QC labs.</p> <p>This includes, but is not limited to</p> <ul style="list-style-type: none"> <li>• Convenient arrangement of functions into four context areas: Sample Entry, Instrument Status, Results Review and Method Management</li> <li>• Graphical instrument status display</li> <li>• "Explorer"-like Tree-View to conveniently display search results from the database</li> <li>• Menus (pull-down as well as context menus)</li> <li>• Toolbars</li> <li>• Tables configurable to the method</li> </ul>
User specific storage of user interface configuration details (profiles)	<p>Specific to job roles and analytical methods.</p> <p>If a user has no permission to use a certain system function, the function is not shown in the user interface or it is disabled (grayed out).</p>

## Analytical instruments

Supported analytical techniques	LC, GC, A/D converter (general purpose)
Number of instruments that system can simultaneously control and acquire data from.	No hard-coded limits (function of database server configuration and available network bandwidth). Typical configurations have 30-50 users and 80-100 instruments. Factory default configurations of Cerity C/S have been setup, optimized and tested for configurations with up to 100 concurrent users and up to 200 instruments. Larger configurations should be planned and optimized based on usage assessment by an Agilent Technologies Laboratory Informatics specialist.

## Data acquisition interfaces

LAN	Level 4 instrument control of 1100 LC, 6890/6850 GC and 35900E ADC, using TCP/IP protocol as installed with Microsoft Windows.
IEEE-488 (HP-IB/GP-IB)	Based on standard interface control library (SICL).
Level of bi-directional instrument control	Level 4 instrument control (advanced full control) for instruments with appropriate capabilities (e.g. Agilent 1100 HPLC, Agilent 1090 HPLC, Agilent 5890 and Agilent 6850/6890 GC). Level 3 instrument control (full control) for Waters 2690/2695 Alliance LC and the Waters 2487 Dual Wavelength Detector. Level 3 (full control) for Shimadzu LC10Avp using the Cerity Generic Instrument Control (GIC) framework.
Instruments for which bi-directional control	Agilent 1100 HPLC, Agilent 6850/6890 GC, Agilent 1090 HPLC Series II (all pump configurations and diode array detector), Agilent 5890 GC (through Cerity Generic Instrument Control framework), Agilent 35900E ADC for general purpose / multi-vendor interfacing., Shimadzu LC10Avp HPLC Waters Alliance 2690 and 2695, Waters 2497 Dual Wavelength Detector is supported.
Proprietary control hardware required?	No. Instruments interfaced via LAN require appropriate instrument LAN interface and an Ethernet LAN card in the PC. Waters Alliance instruments are interfaced via GPIB (IEEE-488) and require an Agilent 82350 GPIB board in the PC.
Network monitoring capabilities (computers and instruments)	The networked design of all components ensures that all parts of the system can be monitored as nodes on the network.

Commercially available network monitoring tools allow to measure network bandwidth, network health, track errors and alerts and help in troubleshooting problems related to communication between computers and instruments on the network.

Network monitoring tools cannot be used to monitor non-networked instruments using legacy connections such as GPIB or RS232.

## Data Transfer

Data import formats	<ul style="list-style-type: none"> <li>• Data files from Agilent ChemStation A.03.01 or greater</li> <li>• ANDI (Analytical Data Interchange) format</li> <li>• Import of work-lists in XML format, e.g. from a LIMS.</li> </ul>
Data export formats	<ul style="list-style-type: none"> <li>• Microsoft Excel format (XLS) for tabular data</li> <li>• HTML for analytical reports</li> <li>• JPEG, GIF, TIFF, WMF for graphics</li> <li>• ANDI (Analytical Data Interchange) format</li> <li>• XML for export of archive catalog files (for use with knowledge management or archive management system)</li> </ul>
Standard interface protocols supported between network components of the system	COM+/DCOM
Results review (o-line on-demand viewing)	By design, all result records are available for instantaneous online review in the results review context of the application.

continued ...



Direct FAX or PDF output	Through standard operating system functions (print to fax, print to PDF device).
Other interfaces	<p>Examples for further customization based on the following techniques</p> <ul style="list-style-type: none"> <li>• Reports post-processing through DOM (Document Object Model)</li> <li>• Reports post-processing through embedded scripts (Javascript or VBScript code)</li> <li>• Web-Access to queries through SSL (Secure Sockets Layer)</li> <li>• Programming interface (requires PSO consulting)</li> </ul>

## Documentation

Description of documentation delivered with the Cerity system	<ul style="list-style-type: none"> <li>• How-To Tasks (Online Help),</li> <li>• Cerity Quick Reference Card</li> <li>• Cerity Concepts Guide</li> <li>• Getting Started Cards</li> <li>• Cerity Installation Guide</li> <li>• Generic Instrument Control</li> <li>• System Administration Guide (Online, PDF)</li> <li>• Cerity Technical Reference Guide</li> </ul>
Description of formats used for online documentation.	PDF and HTML Help
Number of printed manual sets provided with system	<p>One set per user license</p> <p>(Note: Additionally, PDF-versions of all manuals are included on the software media)</p>
Link to standard operating procedures SOPs	The application user interface allows to configure a link directly to the intranet location where SOPs or monograph are stored.

## Access Security and Control

Security concept	<p>Based on NT security system (user accounts management, password policy).</p> <p>Application does not have a proprietary account system but allows reusing password and security policies directly from the operating system.</p>
Access controls for security configuration	<p>User with Cerity system administrator permissions.</p> <p>Note: Cerity system administrators require system administration permissions on the local computer.</p>
Granularity of security access	<p>Managed at the individual function level of the application.</p> <p>Menus and toolbar functions can be selectively configured for each user role. Examples: Create method, reprocess chromatogram, approve result etc.</p>
Description of application security controls	<ul style="list-style-type: none"> <li>• Cerity users must be authenticated through the operating system.</li> <li>• Mandatory login using user-ID and password.</li> <li>• Cerity uses a role-based security concept based on job roles and job responsibilities of users.</li> <li>• Prior to executing a function in the system, Cerity's security service checks whether the user has the appropriate capabilities.</li> <li>• Mandatory audit trail every time a record is created, modified or destroyed.</li> <li>• Cerity allows to configure which system tasks require authorization by electronic signature (e.g. accept/reject an analysis result).</li> <li>• System-wide inactivity time-out locks the session after a predefined idle period.</li> <li>• Physical access security controls are not enforced by the application, but the system is compatible with physical access security controls such as bio-metric or smart-card identification as supported by the operating system.</li> </ul>
Security mechanism of network data packets	Yes, using COM+ security.
Biometrics-based identification	<p>Planned for subsequent release:</p> <p>Support biometrics-based identification with standard interfaces for face recognition, voice recognition, fingerprint scanning.</p>

## Instrument Control, Instrument Status Monitoring and Analysis Scheduling

### Instrument status monitoring

- Graphical display of current online instrument status.
- Transparent real time access to any instrument connected to the network, independent of the instrument and the client computer.
- Transparent access to instrument control and equilibration functions such as reset injector, lamp on/off, set detector wavelength, balance detector, wavelength calibration.
- Real Time Plot.
- Instrument actuals (run-time, instrument errors, warnings, instrument diagnostics data)
- EMF limit checking and reporting for the Agilent 1100.

### Real-time display

Configurable online plot for detector channels and diagnostic plots (e.g. thermostat temperature, pressure, flow).

### Instrument actuals

Configurable status information table (run-time, instrument errors, warnings, diagnosis buffers).

### Real-time status – finding out why the instrument is in an error state or not ready.

Cerity offers three troubleshooting mechanisms for problems such as:

- instrument status GUI and instrument actuals
- instrument logbook
- "Service Report" function that queries the diagnostic registers of the Agilent 1100 to generate a service report that helps to diagnose instrument problems down to the module level. The diagnostic report includes diagnostic information such as instrument error messages (e.g. "leak detected"), firmware revision, EMF limits, pump ripple.

### Scheduling of analyses ("chaining of analyses")

The system uses a scheduling process that allows submitting analysis jobs (single sample analyses as well as sequences) to the so-called worklist. The analysis priority is entered at sample entry and can be high, medium or low.

### Early Maintenance Feedback (EMF)

EMF is an important measure for trustworthy and reliable instrument data in terms of level 4 instrument control. Cerity supports the EMF functions available in the Agilent 1100 HPLC. The EMF concept generates EMF warnings when an user-defined wear and tear limit is reached for one or several parts in the chromatographic system. Cerity logs EMF warnings in the sample logbook. This information can be used to prove that the instrument was in good working order when the original analysis was performed.

### Shutdown methods and standby conditions

Cerity supports so-called "error methods" on the Agilent 1100 HPLC. The instrument will execute the error method after a defined time interval.


In addition, Cerity allows for a so-called "soft-timeout" after which the module will be turned off.

The time interval for initialization of the shut down commands is configurable from 0-999 min.


This approach permits running the instrument in standby-conditions, such as a minimal flow rate with the lamp turned off, and avoid that the system will run out of solvent while being unattended.

In case of a serious instrument error, such as a solvent leak, the system will shut down regardless of the timeout settings to prevent system damage or safety hazards.

A



B



**Figure 3**  
Transparent access to connected instruments, including early maintenance feedback (EMF) information.

## Sample/Sequence Entry

Automatic data entry	A "Sequence Template" can be configured as part of the Certity method. When creating a new sequence, the new sequence is pre-filled based on the settings defined in the sequence template stored in the analysis method.
Mechanisms to minimize typing effort during data entry	<ul style="list-style-type: none"> <li>• Fill-down column</li> <li>• Intelligent fill-down wizard</li> <li>• Apply changes to a selection of sequence lines</li> </ul>
Multiple stop times in sequence methods	The sequence template allows setting analysis times ("stop times") for the different runs of a sequence. This is useful for cases where the analysis time for samples is significantly longer than for calibration standards.
Revision control of sample data	Sample data is subject to strict revision control within the application.
Entry of calibration standards	Weights (concentrations) of calibration standards can be directly entered during sample entry. Partial calibration and multi-vial standards are supported (standard compounds of a certain level can be provided in different vials).
Method specific sample variables	Method specific sample variables (multipliers, divisors) are entered during sample entry and available for custom calculations.
Custom naming of method specific sample variables	Names for sample entry fields are configurable per method.
Naming conventions allow for long descriptors	Names and description fields for samples, sequences, methods and instruments permit long descriptive names with at least 128 characters.
Operator can selectively run data acquisition, processing, reporting	When scheduling an analysis, operators can selectively define which processing tasks will be performed. Options include data acquisition, data analysis and reporting.
Scheduling of analysis	<ul style="list-style-type: none"> <li>• The system uses a scheduling process that allows submitting multiple analysis jobs (single sample analyses or sequences) to the instrument's worklist.</li> <li>• The analysis priority is entered at sample entry. The choices are high, medium or low.</li> <li>• The system permits pausing and resuming of sequences and editing the currently running sequence.</li> </ul>
Intelligent sequencing based on dynamic limit evaluation	<ul style="list-style-type: none"> <li>• The system allows to pause, stop or continue (with warnings) sequences based on automatic evaluation of calculated results against user-defined limits.</li> </ul>

## Method Management

Change control to methods maintained in the data system	Cerity subjects all methods to strict revision control and audit trail logging. No part of the method can be overwritten.
Who needs to modify methods	Method management is typically only required for a senior chemist or chemist. Analysts (chemists, technicians) typically work with predefined methods.
Master methods (template methods)	Master methods are protected methods that serve as template for instrument specific methods. By using a master method, analysts can be sure that the analysis settings are coherent with the official analysis procedure.
Linking methods to analysis procedures	Cerity allows translating a chromatographic analysis procedure into a Certity method. This includes the required sample entry variables, calculations, recalibration schemes, layout and reporting requirements and templates for the sequence of injections.

continued ...


System checks to ensure consistency between method parameters and physical instrument capabilities	<p>The software queries the current configuration details from the physical instrument prior to starting the analysis to ensure that the analysis parameters are compatible.</p> <p>Cerity methods are specific to an instrument configuration (e.g. Agilent 1100 Series LC with DAD and a binary pump) and can be applied for groups of instruments that have the same configuration.</p>
Typical interaction of an analyst with a method	<ul style="list-style-type: none"> <li>• Login to system</li> <li>• Select sample entry context (view)</li> <li>• Login the sample</li> <li>• Select the method (the system automatically suggests the instrument suitable for this analysis)</li> <li>• Enter sample information (description, name, vial number, product code, LIMS ID, concentration of standard etc.)</li> <li>• Schedule analysis</li> </ul>
Settings controlled	<ul style="list-style-type: none"> <li>• Sample variables (numeric input values determined during sample preparation and required for calculations)</li> <li>• Example chromatogram (i.e. a typical chromatogram as generated with this analysis) for display in online results review and on reports</li> <li>• Instrument control/data acquisition</li> <li>• Integration</li> <li>• Peak identification (based on retention time)</li> <li>• Spectral compound confirmation (based on spectral comparison)</li> <li>• Spectral compound purity (based on similarity calculations of compound spectra)</li> <li>• Calibration</li> <li>• Quantification</li> <li>• Custom calculations</li> <li>• Reporting</li> <li>• Data review layout</li> <li>• Reporting Limits</li> </ul>
Reintegration	<ul style="list-style-type: none"> <li>• System allows reintegrating results in the controlled environment of Cerity results review context.</li> <li>• Fine-tuning of integration settings on a specific chromatogram is performed under strict revision control of the result record.</li> <li>• Fine-tuning of integration settings on a specific chromatogram remain private to the chromatogram and do not implicitly affect the master method.</li> </ul>
Reprocessing	<ul style="list-style-type: none"> <li>• Reprocessing functions allow to reprocess data with modified sample parameters, or a different revision of the same method or a different method.</li> <li>• Reprocessing calculations are subject to audit trail and revision control functions of the software.</li> </ul>
Sequence template	<ul style="list-style-type: none"> <li>• The method allows storing the "Sequence Template".</li> <li>• The sequence template defines the normal sequence of injections required for analyses run with this method: blank injections, system suitability, standards, samples, QC samples etc.</li> <li>• The sequence template of the method minimizes data entry effort during sequence setup.</li> </ul>

## Data Analysis

Integration algorithm	Revised version of the Agilent Enhanced Integrator.
Integration events	<p>The system allows setting integration events to change integration parameters appropriate for the signal measured during the analysis. Typical integration events include, but are not limited to:</p> <ul style="list-style-type: none"> <li>• Area reject</li> <li>• Height reject</li> <li>• Slope sensitivity</li> <li>• Peak width</li> <li>• Shoulder detection</li> <li>• Tangent skimming</li> <li>• Detection of negative peaks</li> </ul>
Standard quantification modes	<p>Area%</p> <p>Norm%</p> <p>External Standard Quantification (ESTD)</p> <p>Internal Standard Quantification (ISTD)</p>
Description of recalibration schemes for sequence analyses	<p>The system supports flexible calibration schemes:</p> <ul style="list-style-type: none"> <li>• Moving average calibration (single update calibration)</li> <li>• Standard bracketing</li> <li>• Overall bracketing (also known as "grand average bracketing")</li> </ul>
Description of overall bracketing calibration scheme	<p>Two bracketing modes are available. Overall bracketing calculates one calibration curve per calibrated compound for the sequence and uses it for the quantification of all samples in the sequence. In terms of validation and traceability, this is a lot easier to handle than other floating average recalibration schemes.</p>
Description how the system prevents discrepancies between printed reports and results displayed on the screen.	<p>The Certify report writer is a rendering device that only displays data already stored in the central data repository and does not perform any calculations of its own. This is to avoid discrepancies between results shown on screen and paper.</p>
Description of how system controls re-integration and reprocessing in a controlled manner according to GMP and 21 CFR Part 11.	<p>Every modification of a test result (e.g. in the course of reintegration or reprocessing using updated calibration information) results in a new revision of the result record along with tight links to the metadata and result.</p> <p>Manual intervention in the integration of a chromatogram is notified in the report and in the system's online results view using the following measures:</p> <ul style="list-style-type: none"> <li>• Definition of user capabilities for authorized access</li> <li>• Strict auditing with mandatory audit comments</li> <li>• Authorization by electronic signature (configurable)</li> <li>• Data review layout and Report templates set up to include analysis audit trail and method version</li> </ul>
Support of method specific calculations	<p>The system supports method specific calculations through a spreadsheet control integrated into the data analysis method. The control enables for calculation formulae in Excel-compatible syntax.</p>
Supported calculation functions	<p>Arithmetic, logical, statistical functions also available in off-the-shelf spreadsheet programs.</p>
Storage of method specific calculation results	<p>Yes, including calculation formulae.</p>
Triggering of warnings (pass/fail information) based on calculation results	<p>Yes, such as warnings if system suitability or other limits are exceeded.</p>

Description of user-defined (custom) calculations	<p>The built-in custom calculator is an Excel-compatible spreadsheet control embedded into the Certity data analysis method. It allows setting up method specific calculations.</p> <ul style="list-style-type: none"> <li>• Enables method specific calculations.</li> <li>• Calculations and calculation results are subject to strict revision control within the Certity method.</li> <li>• Eliminate manual transcription and manual recalculation of analysis results.</li> <li>• Custom calculation formulae are entered in Excel syntax.</li> <li>• Custom calculations result formatting conforms to pharmacopoeial requirements when checking results against specifications.</li> <li>• Existing custom calculations can be edited, moved, extended if calculation requirements change (e.g. due to an SOP change).</li> <li>• The Certity custom calculator enables definition of custom calculations on custom calculated results ("calculations on calculations").</li> <li>• Calculation formulae are printable for review and verification.</li> <li>• All custom calculated results are available for online results review and reporting.</li> </ul> <p>Examples for method specific calculations enabled with the Certity custom calculator:</p> <ul style="list-style-type: none"> <li>• Impurity calculations</li> <li>• Reproducibility calculations for replicates</li> <li>• Group statistics</li> <li>• Response factor statistics (calibration precision)</li> <li>• Content uniformity</li> </ul>
Description of system suitability criteria and limits.	<ul style="list-style-type: none"> <li>• Certity calculates peak performance parameters according to the different pharmacopoeias (USP, EP, BP, JP and DAB). The user can configure which of them are reported and shown.</li> <li>• System suitability limits can be defined by component and sample type.</li> <li>• Available noise calculations: peak to peak and ASTM.</li> </ul>
Description of peak identification mechanisms.	<ul style="list-style-type: none"> <li>• The system supports identification by absolute and relative retention times (RRT).</li> <li>• Peak windows for peak recognition are customizable per peak using absolute or relative peak retention time windows.</li> <li>• Peak summing and peak grouping.</li> <li>• The system compensates for retention time variability during analyses using reference peaks for the RT update.</li> <li>• Peak naming is flexible and allows for long names.</li> <li>• Peak confirmation based on spectral data comparison against a prerecorded reference.</li> </ul>
Calibration capabilities	<p>The system allows multilevel calibration with an unlimited number of levels, fixed amount, variable amount and manual response factors.</p>
Reprocessing	<ul style="list-style-type: none"> <li>• All chromatographic runs in a sequence can be reprocessed automatically ("batch-wise").</li> <li>• The system allows to reprocess analysis results after changes have been made to sample information, integration setpoints, calibration information or other setpoint changes.</li> <li>• Existing data may be reprocessed with a newer revision of the original method or with a different method.</li> <li>• Reprocessing efficiency is high because the system will only recalculate modified results; Results not affected by a change are taken from their previous revision.</li> </ul>





**Figure 4**  
Setting up a method-specific statistical group calculation in the "Custom Calculator".

Available calibration schemes

Cerity supports the following calibration curves types: piecewise, linear, quadratic, cubic, exponential, logarithmic, power, average slope.

The following weightings are supported:

- Equal
- # of calibrations
- Linear (x) – by the factor 1/Amount
- Quadratic (x) – by the factor 1/Amount<sup>2</sup>
- Linear (y) – by the factor 1/Response
- Quadratic (y) – by the factor 1/Response<sup>2</sup>
- Lg (x) – by the factor 1/lg(Amount)
- Lg (y) – by the factor 1/lg(Response)
- Ln(x) – by the factor 1/ln(Amount)
- Ln (y) – by the factor 1/ln(Response)

Calibration curve origin treatment

- Include
- Force
- Ignore
- Piecewise (connect)

Calibration review capabilities


- An authorized user can reject individual calibration points manually from a calibration curve; This action is subject to audit trail and requires authorization by electronic signature if configured so.
- Calibration data display includes regression curve, correlation coefficients, confidence intervals (configurable) and relative residuals.

Quantification capabilities


- Response factors are calculated from the calibration result and stored automatically.
- Response factors can also be entered manually per peak.
- Response factors can be updated automatically after performing a re-calibration.
- Results can be calculated using a factor per peak in its calculations.
- Results can be calculated using the same factor for all unknown impurities in its calculations.
- Specified impurities can be calculated using a specific factor.

Automated limit checking

- During sequence execution, the system can be configured to react on limit checks: Pause, abort or continue the sequence with a warning.
- Limits can be defined for single injections and groups of injections, per sample type and per identified peak.
- Limits can be evaluated for all result data, including custom calculated results.
- Multiple limits can be defined for each item to allow for range checking.



**Figure 4.1**  
Spectral compound purity results review. Purity results are calculated automatically and annotated with clear pass/fail results. If the purity check fails the defined acceptance limits, the peak is not quantified.



**Figure 4.2**  
Spectral compound confirmation results review. Confirmation results are calculated automatically and annotated with clear pass/fail results. If a confirmation result is below the defined acceptance limit, the compound is not identified and not quantified.

## Spectral processing

### Interactive spectra handling

- Configurable spectra display
- Defining Wavelength range, background correction, number of spectra and noise range for setup and display of spectral data in results display for UV Spectra handling
- Quick overview of purity and confirmation results with color coded results annotation for each compound and allows configuration of the spectra display in the results context.
- Interactive selection of spectra, reference spectra and spectra ranges
- Display of selected spectra, spectra overlays, spectra differences (residual spectra) and background spectra

### Spectral compound purity

UV Spectral Compound purity: UV purity results are presented similar to confirmation results window displaying both graphical and numerical results. Purity values display with colour coding according to the pass/warning/rejection level of the calculation result. Purity results are available for interactive recalculation for users with access to interactive purity evaluation

### Spectral compound confirmation

UV Spectral Compound confirmation: In results context, compound confirmation information is presented read-only. Users get a graphical result representation of spectra and similarity values plus all numerical results in a table. Confirmation values display with colour coding according to the pass/warning/rejection level of the calculation result. Users can only modify the result display for original, background and overlaid spectra.

## Online Results Review and Approval

Display of analysis results

- Online results review is a separate context (view) of the application..
- Results can be queried from the Cerity database using standard or customize queries (database searches).
- Results are displayed according to the settings defined in the data review layout of the method.

Results approval

- The system supports a 3-step results approval process (analyst review, peer review, manager approval).
- Data can be approved, rejected, or marked for rework.
- Operational system checks ensure the approval steps are performed in permitted sequence of steps.
- Approved data can be locked by the responsible analyst
- Locked data can be unlocked by an authorized administrator

Sample Result	Review Status	Analyst Review	Peer Review	Final Review	Injection Date	Injection Time
1 standard #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:28:14 PM
2 sample1 #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:34:27 PM
3 sample2 #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:40:39 PM
4 sample3 #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:46:51 PM
5 sample4 #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:53:04 PM
6 sample5 #1	Not Done	Not Done	Not Done	Not Done	11/4/2003	5:59:19 PM


**Figure 5.1**  
Result approval in the Accept/Reject screen for results sign-off.

Description of retrieval capabilities

- Database queries can be defined with a query wizard to search for data.
- Queries can be defined based on Samples, Methods, Instruments and results.
- All results can be stored and retrieved.

Protection of electronic records managed by the system

All binary raw data is handled by the Oracle database and Cerity information manager objects and is under strict revision control of the Cerity security service.



**Figure 5.2**  
Cerity results review displaying the results in an Explorer-like treeview. The most recent revision of a record is always displayed first. Historical revisions of the record, including raw data, meta data and audit trail can be viewed by clicking on the "old revisions" mode in the treeview.

The actual position of analyzed vials in auto-samplers is stored and reported (Part 11).

Yes, for Agilent 1100 Series LC and Waters Alliance LC. For 3rd party autosamplers, requires BCD vial number input to the 35900E.

All components of the system are identified in the system (Part 11).

Automatic tracking and storage of instrument serial numbers and firmware revisions (depends on instrument capabilities).


## Reporting

Description of Cerity report generator

The reporting functions consist of a graphical report template editor and a reporting tool that extracts the data to be reported from the Cerity database. Report design is done through drag & drop editing of the report. Report templates are stored in HTML format. Reports can be published to a web browser and on paper.

Configurability of the report generator

- All Cerity reports can be customized using the report template editor.
- Report template editor allows to render (visualize) all data from the application database.
- Item selection is done using a tree view representation of the system's "data dictionary".



**Figure 6**  
Example report template.

Description of built in reports

Single injection reports, group summary, composite (sequence summary) reports. Reports include headers, footers, pagination, logos, result tables, graphics and statistics.

Description of reporting capabilities

- WYSIWYG (report preview in Internet Explorer).
- The report template editor is HTML based.
- Full control of format attributes (layout, fonts, colors).
- Automatic pagination.
- Headers and footers.
- Reports can be printed in 'Portrait' and 'Landscape'.
- Chromatograms can be printed separately. (full page).
- Result files can be printed separately from chromatograms..
- Chromatograms and results can be combined on the same page.
- The system can produce summary tables per analyte.
- Result reports include integration codes per peak.
- Peak markers are shown on the screen and in the report.
- Peak names and retention times can be added to the chromatogram.
- Date and time of printout is reported.
- Date and time of processing or reprocessing is reported.
- Date and time of acquisition is reported.
- The scaling of a chromatogram in a report is user definable.

Conditional reporting (formatting based on limit checks)


The Cerity report template editor allows automatic formatting of report sections such as tables, table columns or individual fields based on logical conditions, e.g. to print a certain result in a specific color depending on its value.

Report printing (after reprocessing)

Reports are automatically generated when reprocessing an analysis e.g. after modifications of sample information, integration settings or other changes to the setpoints.

Report printing (without reprocessing)

- For analysis results already stored in the Cerity database, reports can be regenerated at any time from the results review context without requiring reprocessing.



**Figure 7**  
Example report published in a web browser.

## Data Archiving

Description of Cerity data archiving capabilities

- The built-in archive/restore utility can be used to exchange electronic records (accurate and complete copies) as well as to restore and replay data throughout the record retention period.
- Easy transfer of electronic records to other disks or media for long-term storage and to free up database space.
- Complete audit-trail of all archiving and delete operations.
- Data selection is performed using an archive query wizard.
- XML-based archive catalog allows for interface to archive management tools.

Description of measures to ensure data integrity of archived data

- In order to maintain data integrity, Cerity archives related records in one consistent archive.
- By design, the system prevents archiving incomplete information (e.g. an individual injection from sequence, or injection results without audit trail or method information).
- If archived data was deleted from online storage, it needs to be reloaded to be accessed by the system.

Interface to Information Management Systems

Cerity follows an open system paradigm and can be interfaced to the following types of information management systems:


- Cerity can be interfaced to Agilent Cerity Enterprise Content Management Systems (ECMS) are ideal to manage results data and database archives generated with Cerity and can use Cerity's XML-based archive catalog for keyword indexing.  
Agilent Technologies is a technology partner of Nugenesis Technologies. Agilent Cerity can be interfaced to Nugenesis SDMS.
- Laboratory Information Management Systems (LIMS)  
Cerity provides a generic LIMS interface for the download of LIMS worklists into Cerity and upload of results data into LIMS.  
Agilent Technologies is a technology partner of Labtronics, Inc. Labtronics LimsLink<sup>EDS</sup> provides configurable interfaces for integration with the major chromatography data systems, including Agilent Cerity.
- Enterprise Resource Planning Systems (ERP)  
Cerity for Pharmaceutical QA/QC integrates with BayCovin Instruments. This software from Bayer Business Services bridges the gap between mySAP PLM QM and the laboratory instrumentation level. Using an SAP-certified interface to QM-IDI, BayCovin Instruments is designed to provide a work list of inspections to a client depending on analytical method and instrumentation requirements. Raw quality inspection data is collected from a wide range of instruments, managed, and archived in compliance with FDA regulations (21 CFR Part 11) and ISO quality standards. After optional validation, measured and calculated end results are sent back to mySAP PLM QM to trigger subsequent actions in the logistic chain.

Support of multi-tiered archiving technology

Requires 3rd party tool e.g., active data on hard drive, older data on slower optical magneto-optical disks, jukeboxes.

continued ...

Where are electronic signatures used in the system? What information is stored by the system for each use?	Configurable, depending on the organization's workflow. Includes, but is not limited to, method changes, adding/changing user permissions, result approval.
Are experiments and/or reports reviewed and approved in the system?	Yes. The information can be retrieved and inspected in electronic form as well as on paper for investigations or regulatory inspections.
What features does the system provide to administer user accounts?	Cerity users must be authenticated Windows users. The Cerity system administration console is based on the Microsoft Management Console (MMC) and is used to set the Cerity specific permissions.
Can the system support password aging?	Yes (reuses password policies defined on the operating system level).
Can the system support disabling / re-enabling user accounts?	Yes, by directly reusing the account policy defined in the Windows operating system.
Can the system support account lockouts after a defined number of failed attempts to log-in?	Yes, see above.
How are failed login attempts recorded by the system? How does the system administrator gain access to information about failed login attempts? How are other security problems identified, recorded and accessed by the system administrator?	Standard Windows security, event log and password policies.



**Figure 8**  
Audit trail and electronic signature for results sign-off.




## Functional Specifications – Electronic Records and Electronic Signatures Checklist (21 CFR Part 11)

### Procedures and Controls for Closed Systems

Question	Y/N/NA*	Comments
11.10(a) Is the system validated?	Y	<p>The software development is following the software life-cycle and quality management system of Agilent Technologies Lifescience and Chemical Analysis Group. Full installation qualification (IQ) and operational qualification (OQ) services available for the software as well as the instrumentation controlled by the system.</p> <p>The product was designed to fulfill the validation requirements of the users of this product according to current regulations and quality standards including, but not limited to, 21 CFR 210 (Good Manufacturing Practice for Drugs), 21 CFR 211 (current Good Manufacturing Practice for finished pharmaceuticals), 21 CFR 58 (Good Laboratory Practice), 21 CFR Part 11 (Electronic Records and Signatures).</p>
11.10(a) Is it possible to discern invalid or altered records?	Y	<p>The data managed by the system is under strict access control, revisioning and automatic audit trail. The implementation ensures that alterations to records by authorized individuals result in a new revision of the respective record along with a detailed audit trail. Software operational qualification protocols and database integrity check utilities allow detecting corrupted data (e.g. due to a technical failure of a computer hardware component). Certity generates audit trail entries when a manual integration result is created interactively and discarded later on.</p>
11.10(b) Is the system capable of producing accurate and complete copies of electronic records on paper?	Y	<p>Accurate and complete copies of electronic records are created and handled by the Certity archive/restore utility included in the standard product.</p> <p>In addition, the Certity system provides an HTML based report editor that has the capability to render (visualize) the complete data records managed by the system.</p>
11.10(b) Is the system capable of producing accurate and complete copies of records in electronic form for inspection, review, and copying by the FDA?	Y	<p>Data can be reviewed online in the query-based "Sample", "Instrument", "Method" and "Results Review" views. Navigation through hierarchical data is done using a treeview (similar to Windows Explorer). The Certity archive/restore utility is designed to create accurate and complete copies of electronic records by maintaining complete referential integrity within the archive.</p>
11.10(c) Are the records readily retrievable throughout their retention period?	Y	<p>The system manages all records under strict protection and revision control. Data is accessible directly from the user interface through predefined and customer definable selection criteria (queries). The system implements technical controls so data must be archived before an authorized system administrator can remove it from the online database.</p>
11.10(d) Is system access limited to authorized individuals?	Y	<p>The system requires mandatory login. The system administrator grants user rights (capabilities) to individual users and groups of users. The system checks the user's capabilities prior to executing each task. The system's security implementation uses a combination of operating system (Windows 2000/Windows XP) security/password policies and application specific access controls.</p>
11.10(e) Is there a secure, computer generated, time stamped audit trail that records the date and time of operator entries and actions that create, modify, or delete electronic records?	Y	<p>The system keeps a detailed, human-readable audit trail ("logbook"), which is automatically maintained independently from operators. The audit trail documents every time a record is created, modified or destroyed.</p> <p>Options for electronic sign-off on system tasks, configurable and mandatory audit comments as well as system checks for electronic results review and approval can be enabled using the optional Certity GMP module (Agilent P/N G4030AA).</p>

\*Y=Yes, N=No, NA=Not Applicable



	<b>Question</b>	<b>Y/N/NA*</b>	<b>Comments</b>
11.10(e)	Upon making a change to an electronic record, is previously recorded information still available (i.e., not obscured by the change)?	Y	The system implements strict revision control of its records. No previous entry is ever overwritten.
11.10(e)	Is an electronic record's audit trail retrievable throughout the record's retention period?	Y	The audit trail cannot be separated from the original record. The Cerity archive/restore utility maintains complete referential integrity of the record and archives the record along with its associated metadata and audit trail information.
11.10(e)	Is the audit trail available for review and copying by the FDA?	Y	The Cerity audit trail can be reviewed online (e.g. In results review) and can also be included in printed reports.
11.10(f)	If the sequence of system steps or events is important, is this enforced by the system (e.g., as would be the case in a process control system)?	Y	Operational checks: The system implements technical controls that ensure the permitted sequence of steps. Examples: The system enforces the approval-rejection cycle of results data by the person who generated the data (analyst review), by a 2nd individual (peer review) and final sign-off (final approval). The system's archive-delete cycle requires data to be archived prior to deletion.
11.10(g)	Does the system ensure that only authorized individuals can use the system, electronically sign records, access the operation, or computer system input or output device, alter a record, or perform other operations?	Y	Authority checks: For every transaction, the security services determine whether the currently logged on user has the appropriate authorization based on the access permissions configured for this user or user group
11.10(h)	If it is a requirement of the system that input data or instructions can only come from certain input devices (e.g., terminals) does the system check the validity of the source of any data or instructions received? (Note: This applies where data or instruction can come from more than one device, and therefore the system must verify the integrity of its source, such as a network of weigh scales, or remote, radio controlled terminals).	Y	Device checks: The system performs input verification for data entry fields. Invalid fields are marked red. Agilent instruments automatically detect and record serial numbers and firmware versions. The system stores the hostname of the originating client PC when an electronic record is created or modified. For the 1100 HPLC, the system supports the use of column identifications tags that allow to trace and record analytical column information (e.g. batch number, number of injections, dimensions etc.)
11.10(i)	Is there documented training, including on the job training for system users, developers, IT support staff?	NA	Appropriate user and administrator trainings are available from Agilent Technologies. Customized trainings are available on request. Additional procedural controls are required.
11.10(j)	Is there a written policy that makes individuals fully accountable and responsible for actions initiated under their electronic signatures?	NA	This needs to be addressed by procedural controls in the user's environment.

**\*Y=Yes, N=No, NA=Not Applicable**

	Question	Y/N/NA*	Comments
11.10(k)	Is the distribution of, access to, and use of systems operation and maintenance documentation controlled?	NA	This needs to be addressed by procedural controls in the user's environment.
11.10(k)	Is there a formal change control procedure for system documentation that maintains a time sequenced audit trail for those changes made by the pharmaceutical organization?	NA	This needs to be addressed by procedural controls in the user's environment.

#### Additional Procedures and Controls for Open Systems

	Question	Y/N/NA*	Comments
11.30	Is data encrypted?	N	Not applicable.
11.30	Are digital signatures used?	N	Not applicable.

#### Signed Electronic Records

	Question	Y/N/NA*	Comments
11.50 (a)	Does signed electronic records contain the following related information?	Y	See detailed comments below.
	• The printed name of the signer	Y	System shows the printed name of the signer, date/time (local time and timezone information)
	• The date and time of signing	Y	See previous item.
	• The meaning of the signing (such as approval, review, responsibility)	Y	The meaning of the signature is captured in the context of the function currently executed (e.g. peer review/approval) or through a mandatory comment.
11.50 (b)	Is the above information shown on displayed and printed copies of the electronic record?	Y	Signatures become part of the original record. Changes result in new revisions of records and previous entries are never overwritten. The information is available online and on printed reports.
11.70	Are signatures linked to their respective electronic records to ensure that they cannot be cut, copied, or otherwise transferred by ordinary means for the purpose of falsification?	Y	The signature information becomes part of the original record according to the Cerity database schema. Changes result in new revisions of records and previous entries are never overwritten. In addition, audit trail and e-sig are part of the archived/restored records.

#### Electronic Signatures (General)

	Question	Y/N/NA*	Comments
11.100(a)	Are electronic signatures unique to an individual?	Y	The Cerity security implementation is based on operating system security. This allows to directly reuse the user account system defined by the user's IT operation along with the corresponding security and password policies developed and controlled by the user's organization.
11.100(a)	Are electronic signatures ever reused by, or re-assigned to, anyone else?	Y	The Cerity security implementation is based on operating system security and supports appropriate behavioral controls of the user's organization. For instance, the operating system user account must be disabled but not re-assigned to someone else when the respective individual leaves the organization.
11.100(b)	Is the identity of an individual verified before an electronic signature is allocated?	NA	This must be governed by appropriate company policies.

\*Y=Yes, N=No, NA=Not Applicable

**Electronic Signatures (Non-biometric)**

	Question	Y/N/NA*	Comments
11.200(a)(1)(i)	Is the signature made up of at least two components, such as an identification code and password, or an id card and password?	Y	The Cerity system uses operating system security. Login to the Cerity system and electronic signatures require the user ID and password.
11.200(a)(1)(ii)	When several signings are made during a continuous session, is the password executed at each signing? (Note: both components must be executed at the first signing of a session).	Y	Signing a record always requires entering the user-id and password of that user.
11.200(a)(1)(iii)	If signings are not done in a continuous session, are both components of the electronic signature executed with each signing?	Y	A configurable inactivity timeout prevents impersonation after a defined period without user activity. The currently logged on user must re-enter user ID and password to unblock the system.
11.200(a)(2)	Are non-biometric signatures only used by their genuine owners?	Y	Customer policy has to define, implement and maintain a suitable password policy. Cerity uses Windows security system, allowing reuses of the password policies defined in Windows.
11.200(a)(3)	Would an attempt to falsify an electronic signature require the collaboration of at least two individuals?	Y	Yes, but requires appropriate account and password handling policies in the user's organization and IT environment.

**For Tokens, Cards, and other Devices Bearing or Generating Identification Code or Password Information:**

	Question	Y/N/NA*	Comments
11.300(c)	Is there a loss management procedure to be followed if a device is lost or stolen?	NA	This needs to be addressed by procedural controls in the user's environment.
11.300(c)	Is there a procedure for electronically disabling a device if it is lost, stolen, or potentially compromised?	NA	This needs to be addressed by procedural controls in the user's environment.
11.300(c)	Are there controls over the issuance of temporary and permanent replacements?	NA	This needs to be addressed by procedural controls in the user's environment.
11.300(e)	Is there initial and periodic testing of tokens and cards?	NA	This needs to be addressed by procedural controls in the user's environment.
11.300(e)	Does this testing check that there have been no unauthorized alterations?	NA	This needs to be addressed by procedural controls in the user's environment.

**Electronic Signatures (Biometric)**

	Question	Y/N/NA*	Comments
11.200(b)	Has it been shown that biometric electronic signatures can be used only by their genuine owner?	NA	Not applicable. Cerity for Pharmaceutical QA/QC revision A.02.xx is not delivered with biometric ID devices.

\*Y=Yes, N=No, NA=Not Applicable

## Controls for Identification Codes and Passwords

Question	Y/N/NA*	Comments
11.300(a) Are controls in place to maintain the uniqueness of each combined identification code and password, such that no individual can have the same combination of identification code and password?	Y	Cerity uses the account system, security and password policies defined for the operating system. Therefore, this requires appropriate account and password handling policies in the user's organization and IT environment.
11.300(b) Are procedures in place to ensure that the validity of identification codes is periodically checked?	Y	Yes, see 11.300(a)
11.300(b) Do passwords periodically expire and need to be revised?	Y	Yes, see 11.300(a)
11.300(c) Is there a procedure for electronically disabling an identification code or password if it is potentially compromised or lost?	Y	Yes, see 11.300(a)
11.300(d) Is there a procedure for detecting attempts at unauthorized use and for informing security?	Y	Yes, see 11.300(a). The Cerity security implementation uses the operating system event viewer to log security events. This requires appropriate configuration of the operating system's event logging.
11.300(d) Is there a procedure for reporting repeated or serious attempts at unauthorized use to management?	Y	Yes, see 11.300(a)

\*Y=Yes, N=No, NA=Not Applicable

## Warranty and Support Contracts

Please contact your local support sales representative.

Warranty Period	Varies by country and can be from 1-3 years.
Extended Warranty	Available.
Software telephone support	Available.
Software materials subscription	Available.
Software status bulletins	Available from <a href="http://www.agilent.com/chem/nds">www.agilent.com/chem/nds</a> (requires valid software license number)





Figure 9  
Computer based protocol for OQ/PV on an Agilent 1100 Series HPLC system.

Services	Please contact your local NDS sales representative.
Standard services	Installation, Familiarization, Education
Available standard courses	H2296A – Certity Networked Data System Basic Operation (Full Access Users) 1 day H2297A – Certity Networked Data System Advanced Operation (Full Access Users) 2 days H2298A – Certity Networked Data System Application Administration 1 day H2295A – Certity Networked Data System Routine Operation 1 day Customized education courses are available and can be delivered at a central Agilent location or on-site.
Personalized education services	
Qualification services for NDS software	Computer-based installation qualification (IQ) Computer-based operation qualification (OQ/PV) Delivered by Agilent customer engineers or a certified support provider.
Qualification services for chromatography equipment	Computer-based installation qualification (IQ) Computer-based operation qualification (OQ/PV) Delivered by Agilent customer engineers or certified support provider.
Customization services	Available. Delivered through Agilent's Project Services Organization. Please contact your local NDS sales representative.
Qualification of computer network infrastructure	Regulatory agencies are endorsing a risk-based approach to compliance. Network infrastructure can be a high-risk, high-impact component for the integrity and security of data it transports. Network infrastructure needs to be considered in a company's overall master validation plan.  Agilent Technologies has taken a metrology-based approach to the qualification of computer network infrastructure and provides services the following services: <ul style="list-style-type: none"> <li>• Design Qualification for networks</li> <li>• Installation Qualification for networks</li> <li>• Operational Qualification for networks</li> </ul> These services are delivered through Agilent's Project Services Organization.
Consulting services	Available. Delivered through Agilent's Project Services Organization.
Project management services	Available. Project management according to development lifecycle. Services comprise specification, design, implementation, deployment, validation and support. Delivered through Agilent's Project Services Organization.
Declaration of system validation audit reports	Included with the shipment kit.  The quality management system of Agilent Technologies and the lifecycle documentation for Certity for Pharmaceutical QA/QC has been audited by independent inspectors according to PDA Technical Report #32. The audit report is available to subscribers from the audit repository center (ARC) <a href="http://www.auditcenter.com">www.auditcenter.com</a>



**Figure 10**  
The software qualification protocols are based on the Certity TestManager, an automated regression test utility.



**Figure 11**  
The audit report is available to subscribers from the audit repository center (ARC)  
[www.auditcenter.com](http://www.auditcenter.com).



**[www.agilent.com/chem/cds](http://www.agilent.com/chem/cds)**

The information in this publication is subject to change without notice.

Microsoft® and Microsoft Windows® are U.S. registered trademarks of Microsoft Corp.  
Oracle® is a U.S. registered trademark of Oracle Corporation, Redwood City, California.

Copyright © 2003-2004 Agilent Technologies, Inc.

All Rights Reserved. Reproduction, adaptation or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Published July 1, 2004  
Publication Number 5989-1425EN



**Agilent Technologies**