

# Effect of Molecular Weight on Sample Loading in Gel Permeation Chromatography

## Technical Overview

### Introduction

In gel permeation chromatography and size exclusion chromatography, sample loading is very dependent on the molecular weight distribution of the polymer. Broader distribution polymers can be loaded to a greater extent. This effect is demonstrated on an Agilent PLgel 10  $\mu$ m MIXED-B column. Figure 1 shows the column calibration using Agilent EasiCal polystyrene calibrants. Figure 2 shows the effect of sample loading on the molecular weight distribution of a polystyrene with  $M_w = 250,000$  and dispersity = 2.5.

### Conditions

Calibrants	EasiCal (0.1% solution, 200 $\mu$ L injection)
Column	Agilent PLgel 10 $\mu$ m MIXED-B, 25 $\times$ 300 mm (p/n PL1210-6100)
Eluent	THF
Flow rate	10 mL/min
Detector	RI
System	Agilent PL-GPC 50



**Agilent Technologies**

Mp (upper trace)	3,040,000	Mp (lower trace)	8,500,000
	330,000		1,030,000
	66,000		156,000
	9,200		28,500
	580		3,250




Figure 1. Calibration of an Agilent PLgel 10 µm MIXED-B column with Agilent EasiCal polystyrene standards.




Figure 2. Effect of polystyrene molecular weight on sample loading demonstrated on an Agilent PLgel 10 µm MIXED-B column.

## **GPC/SEC columns and calibrants from Agilent**

Agilent offers a comprehensive portfolio of GPC/SEC columns and calibrants for high-performance separations based on molecular size in solution. Agilent delivers leading solutions for characterizing and separating polymers by GPC/SEC, and manufactures all components for accurate polymer analysis.

Look at the Agilent Literature Library on [www.agilent.com/chem/gpc-sec](http://www.agilent.com/chem/gpc-sec) for a comprehensive range of application notes and technical overviews to help you get the best from your Agilent GPC/SEC columns and instruments.

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2011  
Printed in the USA  
July 12, 2011  
5990-8611EN



**Agilent Technologies**