




Efficient Separation of Polymethylmethacrylate

Agilent PLgel 10 μm MIXED-B Columns

Technical Overview

Introduction

PLgel 10 μm MIXED-B columns are designed for high MW polymer analysis and demanding eluent conditions. The PLgel 10 μm MIXED-B spans a wide range of molecular weights, up to 10 million, with a linear calibration curve. It is particularly useful for molecular weight distributions where slightly higher than average MWs are encountered. The 10 μm particle size provides good resolution with relatively low pressures for enhanced lifetimes in demanding conditions.



Agilent Technologies

The figure shows the separation of a polymethylmethacrylate obtained using PLgel MIXED-B columns. By applying Mark-Houwink parameters (polystyrene $K = 14.1 \times 10^4$ $a = 0.70$; polymethylmethacrylate $K = 12.8 \times 10^4$ $a = 0.69$)¹, the differences in molecular weight averages shown in Table 1 are revealed.

Conditions

Columns: 2 x PLgel 10 μ m MIXED-B,
300 x 7.5 mm (p/n PL1110-6100)

Eluent: Methyl ethyl ketone (MEK)

Flow Rate: 1.0 mL/min

Detection: RI




Figure 1. Separation of a polymethylmethacrylate obtained using PLgel MIXED-B columns

Table 1. *M_p*, *M_w* and *M_n* equivalents for PS and PMMA

	PS Equivalent	PMMA Equivalent
M_p	37521	42286
M_w	36669	41628
M_n	19354	21714

¹ B J Hunt and S R Holding (eds) (1989) Size Exclusion Chromatography. Blackie, UK.

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 23, 2010

SI-01778



Agilent Technologies