

Successful Polyetherimide Analysis Agilent PLgel 5 μ m MIXED-D Columns

Technical Overview

Introduction

Polyetherimide elutes well within the recooling range of PLgel MIXED-D columns. The PLgel 5 μ m MIXED-D column is specifically designed for the analysis of polymers, paints and resin systems where material above 400,000 MW is unlikely to be present. High pore volume, concentrated in this operating range, combined with the 5 μ m efficiency, provides excellent resolution for low MW polymers and oligomers. Two, or even three, PLgel 5 μ m MIXED-D columns are the perfect replacement for the popular $10^4/500$ Å or $10^4/10^3/500$ Å/ 100 Å column combinations.



Agilent Technologies

Figure 1 shows the separation of a polyetherimide. The sample was dissolved by heating to approximately 70 °C; this temperature was also used for the analysis due to the high viscosity of the eluent.

Column: 2 x PLgel 5 µm MIXED-D, 300 x 7.5 mm
(part number PL1110-6504)
Eluent: DMF + 0.1% LiBr
Flow Rate: 1.0 mL/min
Detection: RI




Figure 1. Separation of a polyetherimide using a PLgel MIXED-D column set

www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2010

Published in UK, September 21, 2010

SI-01831



Agilent Technologies