

Brilliant® II QPCR and QRT-PCR Reagents

The Next Generation of QPCR Reagents

- + Improved sensitivity and flexibility
- + Earlier Ct detection at high and low target concentrations
- + Increased SYBR® Green fluorescence with shorter run times

OUR NEXT GENERATION BRILLIANT II® QPCR AND QRT-PCR MASTER MIX KITS OFFER SUPERIOR SENSITIVITY FOR IMPROVED QUANTIFICATION AND REPRODUCIBILITY.

Our new Brilliant® II QPCR and QRT-PCR reagents^{a,b,c} were developed to maintain the high quality of our original Brilliant product line, but with a significant improvement in sensitivity of detection, ensuring reproducible quantification even at low target concentrations. When compared to other QPCR reagents on the market, we observe earlier threshold cycle (Ct) detection at high and low template concentrations as well as improved reaction efficiency across the entire range of template concentrations. Earlier Ct detection and higher efficiency of amplification ensure greater reproducibility within an assay and across multiple assays and different templates.

The Brilliant II QPCR and QRT-PCR reagents have been developed, tested, and validated on our Mx3000P® and Mx3005P® QPCR Systems. By working with both our Mx™ QPCR systems and our new

Brilliant II QPCR reagents, you will benefit from utilizing a fully integrated system optimized to function at peak performance.

Improved Sensitivity of Detection

The new Brilliant II QPCR and QRT-PCR master mix reagents are designed to provide you with improved sensitivity when compared with other commercially available QPCR reagents. In Figure 1, our new Brilliant II SYBR Green QPCR master mix outperforms four other competitor SYBR Green based master mixes. The improved sensitivity using the Brilliant II reagents is evident with earlier Ct detection of 1-3 cycles at each template concentration, improved replicate reproducibility, and higher amplification efficiency across four orders of magnitude.




Figure 1
Brilliant® II SYBR® Green QPCR Master Mix Delivers Better Sensitivity and Reproducibility over a Wide Range of Concentration

Standard curve plot of 10-fold dilutions detecting a GAPDH target (150bp, 52% G/C) from 100 ng down to 0.01 ng of cDNA. The Brilliant® II SYBR® Green QPCR Master Mix detects template at ~1-3 Ct's earlier than competitor reagents and with perfect efficiency and Rsq values.

Brilliant® II QPCR and QRT-PCR Reagents

In Figure 2, we illustrate the improved replicate reproducibility of Brilliant II SYBR Green 1-Step QRT-PCR master mix compared to three competitor reagents. The earlier Ct detection and improved replicate reproducibility leads to more accurate quantification of starting material.

Another measure of sensitivity is the ability to reproducibly distinguish small differences in template concentrations. One of the most challenging assays in QPCR is detecting 2-fold differences in samples at very low copy number (or concentration). Our Brilliant II QPCR reagents are capable of quantifying 2-fold differences (equal to 1 cycle or 1 Ct difference) in samples below 10 copies of template (Figure 3). In this example, we quantify the difference between 2.5 copies and 5 copies of target and illustrate 95% amplification efficiency with consistent detection of 2-fold differences at each concentration.

Brighter, Faster SYBR Green QPCR

Our new Brilliant II SYBR Green QPCR master mix exhibits earlier Ct detection with over twice the fluorescence intensity compared to our original Brilliant SYBR Green product and higher fluorescence than many competitor reagents. We also optimized cycling conditions for more rapid run times of 2 hours total (over a 45 minute reduction in overall run time) with increased performance compared to competitor's recommended longer duration cycling conditions.

More Efficient QRT-PCR

An additional upgrade in our Brilliant II QRT-PCR reagents is the use of new reverse transcription (RT) enzymes specifically formulated for improved performance in one- and two-step formats.

The Brilliant II QRT-PCR one-step reagents use a Moloney-based RT specifically formulated for the Brilliant II 1-step kit with optimal performance at a synthesis temperature of 50°C.

The Brilliant II QRT-PCR two-step reagents include our new AffinityScript™ Multiple Temperature Reverse Transcriptase*. The AffinityScript RT contributes to improved sensitivity, but is also designed for greater flexibility to generate high cDNA yields at synthesis temperatures of 42-55°C.




Figure 2
Improved Reproducibility with Brilliant® II Reagents

Box and Whisker chart illustrating the Ct variation of Brilliant® II reagents and three other commercially available kits. SYBR® Green I dye was used for detection of ARF target (131 bp, 54% G/C) gene across 8 replicates at 20 pg of total RNA. The Brilliant® II SYBR® Green QRT-PCR 1-Step Master Mix has min-max delta Ct of 0.85 cycles compared to delta Ct values of 1.2-3.2 (max delta Ct is nearly a 10x fold variation) for competitor reagents.





Figure 3
Superior Sensitivity and Precision to Detect Target Across 2-Fold Dilutions

Amplification plot (semi-log) and standard curve of 2-fold dilutions of linearized plasmid detected with a linear hydrolysis probe. Brilliant® II QPCR reagents exhibit precise detection of 2-fold differences from 1,280 copies down to ~2.5 copies. The average delta Ct across all 10 concentrations is one cycle and the standard curve efficiency equals 95%.

Ordering Information:

Find an Agilent customer center in your country:
www.agilent.com/chem/contactus

U.S. and Canada:
1-800-227-9770
agilent_inquiries@agilent.com

Asia Pacific:
adinquiry_aplsc@agilent.com

Europe:
info_agilent@agilent.com

For a list of worldwide distributors, please visit our website.

For more information please visit our website:
www.stratagene.com

Brilliant® II QPCR and QRT-PCR Reagents and Related Products

Probe-Based Master Mixes

Brilliant® II QPCR Master Mix 400 rxn 600804

10 x 400 rxn 600815

Brilliant® II QRT-PCR 1-Step Master Mix 400 rxn 600809

10 x 400 rxn 600818

Brilliant® II QRT-PCR AffinityScript™ 2-Step Master Mix 400 rxn 600827

SYBR® Green Master Mixes

Brilliant® II SYBR® Green QPCR Master Mix 400 rxn 600828

10 x 400 rxn 600831

Brilliant® II SYBR® Green QRT-PCR 1-Step Master Mix 400 rxn 600825

10 x 400 rxn 600826

Brilliant® II SYBR® Green QRT-PCR AffinityScript™ 2-Step Master Mix 400 rxn 600834

Probe-Based Core Kits

Brilliant® II QRT-PCR 1-Step Core Kit 400 rxn 600810

10 x 400 rxn 600819

cDNA Synthesis for QPCR

AffinityScript™ QPCR cDNA Synthesis Kit 400 rxn 600559

Brilliant®, Ma3000P®, and Ma3005P® are registered trademarks of Stratagene, an Agilent Technologies Company, in the United States.

AffinityScript® and Mx™ are trademarks of Stratagene, an Agilent Technologies Company, in the United States.

SYBR® is a registered trademark of Molecular Probes, Inc.

a. Practice of the patented 5' Nuclease Process requires a license from Applied Biosystems. The purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research when used with the separate purchase of Licensed Probe. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94044, USA.

b. Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94044, USA.

c. SYBR® Green is licensed for research and development only under patents owned by Invitrogen Corporation.

d. Patents pending.

© Stratagene, an Agilent Technologies Division 2008

Printed in the USA, August 2008

5989-8271EN

STRATAGENE

An Agilent Technologies Division