English			
Disposabile ECG cable			
User Manual	Version: A0		

1. Indication for Use/Intended Use

The disposable ECG (electrocardiograph) cable is intended to be used with ECG monitors & recorders or other physiologic recorders that include an electrocardiographic signal in their function. The cable is used to connect electrodes, placed at appropriate sites on the patient, to the aforementioned devices for general monitoring and/or diagnostic evaluation by health care professional.

2. CAUTION

Law restricts this device to sale by or on the order of a physician.

Item	P01-001-00565
Classification	Lead wire
Length (in/mm)	72/1830
Lead number	5
Distal connector design	Snap

4. Operation Environment

- Environment temperature: 0°C ~40°C
- Relative humidity: 15% ~85%
- Atmospheric pressure: 86kpa ~106kpa

5. Storage

Storage conditions are as follows:

- Temperature: -10 to +40°C
- Relative humidity: 15 to 85%
- Atmospheric pressure: 86kpa ~106kpa

6. Safety & Performance Compliance

Safety: IEC 60601-1: 2012 and

Performance: ANSI AAMI EC53: 2013

Biocompatibility: ISO 10993-5: 2009, ISO 10993-10: 2010

7. Installation for use:

To obtain instructions on the electrode placement, refer to the information indicated by the host device.

Warning: MR-unsafe!

Do not expose the device to a magnetic resonance (MR) environment.

- Thermal injury and burns may occur due to the metal components of the device that can heat during MR scanning.
- The device may not function properly due to the strong magnetic and radiofrequency fields generated by the MR scanner.

8. Replacement and disposal





- Replacement: Check the appearance before use. If the cable/jacket damaged when opening the packing bag, please replace the new ECG cable.
- Disposal of the cables should comply with local regulation.

9. Warning


- This cable is for use only with legal marketed host device.
- Please always remember to read the user manual for the compatible monitor before using the ECG cable.
- Do not allow any liquid to get inside the various connectors on the interface cable.
- Do not reuse the disposable ECG cables.
- Do not use a cable or part of a cable if there is any risk to the patient (e.g. damaged insulating material).
- This disposable product is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the product after patient use. Multi-patient use without disinfection between patients increases the risk of infection. This disposable product is not designed for reuse. Change cannot guarantee the performance specifications once this product has been disinfected or cleaned after patient use.
- If, in relation to the use of the cable, a death or a serious deterioration of health has occurred, this should be reported to the manufacturer and the competent authority of your country.
- This cable [or leadwire or set of leadwires] is not defibrillator-proof and should not be used in areas where a PATIENT might be defibrillated.
- This cable [or leadwire or set of leadwires] meets the requirements of the standard, ECG TRUNK CABLES and PATIENT LEADWIRES (ANSI/ AAMI EC53).

10. Contra-Indications

- No novel features.

	Caution		Manufacturer		Date of manufacture
	Type BF applied part		Reference number		Batch code
	Do not re-use		Consult instructions for use		The device should be sent to the special agencies according to local regulations for separate collection after its useful life.
	EU representative		MR-unsafe		Federal (U.S.) Law restricts this device to sale by or on the order of a physician
	CE marking		Medical device indicator		Biodegradable materials

 Shenzhen Changke Connect Electronics Co., Ltd. Address: A1-2nd and A2-4th floor of factory building, No.87 of Hengping Road, Baoan Community, Yuanshan Street, Longgang District, Shenzhen, China

 Wellkang Ltd. Enterprise Hub, NW Business Complex,1 Beraghmore Rd., Derry, BT48 8SE, Northern Ireland, UK Exclusively marketed by Philips

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