剄 FUJIFILM

Image Reader FLA-7000

Operation Manual



Version 2.0

2006. Mar

2006. Mar

Contents

Part 1 Preparations Before Use1
1. Starting the Image Reader 2 2. About the Main Window 4
Part 2 Reading Imaging Plates5
1. Setting the Reading Conditions 6 2. Setting the IP on the IP Stage 8 3. Setting the IP Stage on the FLA-7000 8 4. Starting Reading 9
Part 3 Reading Fluorescent Samples13
1. Setting the Reading Conditions 14 2. Setting a Fluorescent Sample on the FLUOR Stage 17 3. Setting the FLUOR Stage on the FLA-7000 17 4. Starting Reading 17
Part 4 Reading Digitized Samples21
1. Setting the Reading Conditions222. Setting a Digitized Sample on the FLUOR Stage253. Setting the FLUOR Stage on the FLA-7000254. Starting Reading25
Part 5 Lasers and Filters, Other Settings
1. Registering Laser and Filter Combinations (Method) 30 2. Filter Module Settings (Filter Module) 31 2a. Registering Filters in the Image Reader 31 2b. Saving the Four Types of Filter Combinations 34 2c. Recalling Four Types of Filter Combinations 35 2d. Registering a New Filter Name 35 3. Other Settings (Preference) 36 3a. Settings for Scanning (Scan Settings) 36 3b. Selecting the File Format for Saving and the Analyzing Software to Launch (Image File Settings) 38
Part 6 Installing and Uninstalling the Software
1. Installation (For Windows®) 42 1a. Installation of FUJI USB Control Driver 42 1b. Installation of FUJI USB Function Driver 47 1c. Installation of Image Reader FLA-7000 Software 48 2. Uninstallation (For Windows®) 49

 Installation (For Macintosh[™]) Uninstallation (For Macintosh[™]) 	
Part 7 Troubleshooting	51
1. Errors	
2. Warnings	

Preparations Before Use

1. Starting the Image Reader

Caution

Please note that the GUI screens may change without notice.

Note1:

The Image Reader FLA-7000 software is available in two types: a Windows version, and a Macintosh version. Both versions have the same functions.

This manual uses the screens of the Windows version. If you are using the Macintosh version, follow the instructions in this manual, except those of OS-related operations (such as starting and exiting the software).

Note2:

The following computers are compatible:

OS: Windows XP Professional SP2 or later, Mac OS X 10.3 or later Memory: 512 MB Interface: USB 2.0

Precautions for Use

Do not connect any USB devices other than the FLA-7000 to the computer in which the Image Reader is installed. Otherwise, it may cause malfunction.

During reading, do not use any USB devices other than the FLA-7000 connected to the computer. If USB devices are used simultaneously, image data may be lost.

1-1 Turn on the FLA-7000 and peripheral devices.

Caution

If an Imaging Plate is inserted in the FLA-7000 before turning it on, its scanned data cannot be guaranteed. The sensitivity of the Imaging Plate may deteriorate, based on the self-diagnosis of the FLA-7000.



1-2 Turn on the computer.

1-3 Make sure that the FLA-7000 has completed the warm-up. (Only the power lamp on the upper left panel on the front of the FLA-7000 is lit up when the warm-up is completed.)
 Then, perform the procedures below.

For the Windows version: Start the Image Reader FLA-7000 software from the Startup menu, or use the shortcut key.



For the Macintosh version: Double-click the alias or the software to start the Image Reader FLA-7000 software.



Software



Alias

Note:

After starting the Image Reader, its condition is displayed in the Status area.

- When disconnected: Cannot recognize FLA-7000. Please check connection and power.
- During warm-up:
- FLA-7000 is in self-diagnosis. Please wait.
- When reading is possible:
- Ready
- 1-4 The main window of the Image Reader FLA-7000 software is displayed.



2. About the Main Window



4

Part 2

Reading Imaging Plates

- 1. Setting the Reading Conditions
- 1-1 Click the button.The Reader Settings window for the IP mode is displayed.

		5	6	Ć	D
	Image Reader for FLA-7000				
	<u>File Edit H</u> elp				
	Condition File :		Save	Condition	Load Condition
Ul	-Image Folder : ○ ¥Documents and Settings¥Administrator¥ቻ	スクトップ		•	browse
<u>(2</u>)—					
õг	Comment : dfs				
3	Reader Settings : IP Mode				
4)—	Laser : 💽 LD650 Filter : 🔄 [IP]	Sensitivity S10000 S4000 S1000	I Pixel Size ⊂ 25 µm ⊙ 50 µm ⊂ 100 µm	Lat ເ	ı itude L4 L5
(8)—	Sampling Area		⊂ 200 µm		
	ABCDEFGHIJK 1 2 3 4 5 6 7 8 ABCDEFGHIJJK		Mode	out 0 m 40 s out 15.26 ME	3
	D Top Image File Format : Img/Im Image Data Type : Log ND Filter : Off	f		► S	tart Scan

I Top : Return to the main window from the IP mode.

1-2 Make these settings before reading IPs. Refer to the following explanations of reading conditions when making settings.

1:

Image Folder: C*Documents and Settings*Administrator*デスクトップ
 Specify where to save the file after reading.
 Click the browse... button.

2 File Name : test

Input the name of a file for saving data of a read image. You may not start reading unless you input a file name.

③ Comment : dfs :

The comment is saved with the image as a file, and can be viewed with the analyzing software. Input it if necessary.

④ Before starting reading, make sure that the LD 650 nm laser is shown and the IP filter is selected in the Setting field as shown below.

In IP reading, the LD 650 nm laser and IP filter are selected automatically.



Note:

To use the IP mode, the LD 650nm laser must be loaded, and the IP filter must be set. If these conditions are not fulfilled, the Reader Settings window for the IP mode cannot be accessed.



(6) Pixel Size c 25 μm
 c 25 μm
 c 50 μm
 c 200 μm
 <lic 200 μm
 <lic 2

200	100	50	25 µm
Short	Reading	g Time	Long
Small	Image F	ile Size	Large

The dynamic range of the dynamic range. The dynamic range that can be detected is bigger with L5 than with L4. If the signals of the sample are in the L4 range, the density gradation is represented more finely if L4 is selected.



Drag and select the scanning area on the IP stage.

ſ

	Note:		
	Mode Select this to specify the reading area based on the 2.5 cm grid lines on the IP stage. 20 x 25 Select this to specify the reading area arbitrarily. 20 x 40 Select this to specify a reading area of 20 cm x 25 cm. Select this to specify a reading area of 20 cm x 40 cm.		
	Save Condition : Use this button when saving the reading conditions in a file. Yo may save reading conditions that are used frequently with this function and recall them later with Load Condition		
	Load Condition : Use this button when recalling reading conditions saved with Save Condition		
	Note: When starting the Image Reader, the settings information from the previous session is displayed.		
2. Setting the IP on the IP Stage	Set the IP on the IP stage. For instructions on setting the IP, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.		
3. Setting the IP Stage on the FLA-7000	Set the IP stage on the FLA-7000. For instructions on setting the IP stage on the FLA-7000, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.		

4. Starting Reading



The scanned area is displayed in the real-time window, as shown below.

The stage is read from the left towards the right.



Reading may be finished at any time before the whole scanning area has completed reading. Click the <u>Stop</u> button when you want to finish reading.

Note1:

If you click Stop during reading, the part that has not been read yet will be saved as an image with a data value of 0 (light intensity of 0).

Note2:

When you click the Stop button, reading itself is canceled. You cannot start reading again from the location where reading stopped. 4-2 When you want to change the display parameters of the real-time window, refer to the explanations below.

Curve: ☑ Auto	Exponential 💌 Exponential Linear	: Select the ty down menu. Exponential
	Sigmoid	
		Linear:

Select the type of tone curve from the pulllown menu.

Exponential: The exponential tone curve is used to adjust gradations. Linear: The linear tone curve is used to adjust gradations. Sigmoid: The sigmoid tone curve is used to adjust gradations.



Drag the adjuster.

You may adjust the density of the read image.

Data of a lower light intensity that the line on the left (Low value) will be displayed as a completely white image, and data of a higher light intensity than the line on the right (High value) will be displayed as a completely black image.

Auto Range Scope : If Auto Range Scope is checked, the Image Reader automatically corrects the optimum tone.

Note:

The Image Reader converts data read from samples to images that have an information of 65536 tones, with 0 being the value for white, and 65535 being the value for black.

The tones are indicated by the horizontal axis of the tone curve graph.

Magnification	: 10:1 (10%) 💌
	10:1 (10%)
	5:1 (20%)
	2:1 (40%)

: You may change the display area by selecting a magnification ratio from the pull-down menu. In addition, after reading, if you click the display area after clicking the 🗨 🔍 buttons, the clicked area can be enlarged or reduced.

4-3	Save as : Save the data with a different file name.			
	Launch : Launch the registered analyzing software to display the image.			
4-4	To read a second IP continuously, carry out reading by following the above procedures. Click the <u>Return</u> button to return to the first Reader Settings window. Do not open the stage door of the FLA-7000 until the stage has completely returned. If it is opened, close it immediately. When scanning finishes, the <u>Save as</u> . and <u>Launch</u> buttons become active, but the <u>Return</u> button is grayed out until the stage has completely returned.			
4-5	Finish reading.			

Before turning off the power of the FLA-7000, shut down the Image Reader software.

Note:

If the photo-multipliers (PMT) are replaced with a multi-alkali PMT, the IP mode switches from that of the regular IP mode to a PMT voltage adjustment mode. The PMT voltage adjustment mode differs from the regular IP mode in that a sensitivity level that matches the scanned sample can be set independently. Replacement of the standard PMT with a multi-alkali PMT is done by a serviceman. For details, contact an authorized dealer.

The Sensitivity settings differ from that of the regular IP mode. The operation procedures for other settings are the same as those of the regular IP mode. Follow the operation procedures of the regular IP mode.

IP mode					
Sensitivity	Pixel Si ⊂ 25 ← 50 ⊂ 100 ⊂ 200	ze μm μm μm μm	Latitude ເ⊂L4 ⊂L5		
PMT voltage adjus	tment n	node			
PMT	Pixel Siz	ze	Latitude		
500 V (200-1000)	○ 25	μm	⊙ L4		
	⊙ 50 ○ 100 ○ 200	μm μm μm	⊂ L5		
You may	set the	voltage	to be appli	ed to the	PMT as an
integral v	alue wit	hin the p	oredetermi	ned range	e. The larger
the value	is, the l	higher th	ne sensitivi	ty, but no	ise will be
areater	,	0			
grouton		200	V	Malua	1000V
		Small	1	value	Large

Low <

Sensitivity

High

Part **3**

Reading Fluorescent Samples

- 1. Setting the Reading Conditions
- 1-1 Click the Fluorescence button.
 The Reader Settings window for the Fluorescence mode is displayed.



I Top : Return to the main window from the Fluorescence mode.

- 1-2 Make these settings before reading fluorescent samples. Refer to the following explanations of reading conditions when making settings.
 - Image Folder: C+Documents and Settings #Administrator サデスカトップ
 Click the browse... button and specify where to save the file after reading.

Note:

When scanning multiple times, the image data is saved in the folder with the name specified in File Name, which is automatically created in the specified location.

2 File Name : test1 :

Input the name of a file for saving data of an image.

You may not start reading unless you input a file name.

<i>Note:</i> When the number of scans is between 2 and 4, the laser name and scan number are automatically added to the specified name.
000- 635- 1 File Laser Scan name number
Example: If reading is set with the LD 635 nm laser, LD 473 nm laser, LD 635 nm laser, and LD 635 nm laser, with "test" as the file name, then the following files are created in the "test" folder. test-635-1 test-473 test-635-2 test-635-3

3 Comment : test

The comment is saved with the image as a file, and can be viewed with the analyzing software. Input it if necessary.



3rd : Scanning can be performed up to 4
 times. The following items can be set.

:

Note:

Method : MT :	[Cy5] [SYBR Green] [FITC] [Cy2] [Cy5] [IP] [Digit473]	From the pull-on that correspon selected laser displayed belo	down menu, select ds with the sample and filter combina w.	the Method The tion is
PMT :	500 V (200-1000)	: You may set the within the predete value is, the higher	roltage to be applie be (PMT) as an in rmined range. The er the sensitivity.	ed to the tegral value a larger the
		200V Small ◀ Low ◀	Value Sensitivity	1000V ──► Large ──► High
-	+ : Click th and clic Up to 4	he + button to in ck the - button t 4 scans can be perfo	crease the numbe o reduce the numl ormed.	r of scans, per of scans.

Pixel Size c 25 μm c 50 μm c 100 μm c 200 μm
 Set the pixel size for reading. Click to select from one of the 4 types, as shown on the left. A sample with a smaller pixel size is analyzed more finely.

200	100	50	25 μm
Short	Reading Time		Long
-			
Small	Image Fi	ile Size	Large

(6) Latitude C L4 C L5 : Specify the dynamic range. The dynamic range that can be detected is bigger with L5 than with L4. If the signals of the sample are in the L4 range, the density gradation is represented more finely if L4 is selected.



Size of the file. Expressed as file size per test x number of scans.

Drag and select the scanning area on the FLUOR stage.



Save Condition... : Use this button when saving the reading conditions in a file. You

	may save reading conditions that are used frequently with this function and recall them later with LondCondition
	Load Condition : Use this button when recalling reading conditions saved with Save Condition
	Note: When starting the Image Reader, the settings information from the previous session is displayed.
2. Setting a Fluorescent Sample on the FLUOR Stage	Set the fluorescent sample on the FLUOR stage. For instructions on setting the fluorescent sample, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.
3. Setting the FLUOR Stage on the FLA-7000	Set the FLUOR stage on the FLA-7000. For instructions on setting the FLUOR stage on the FLA-7000, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.

- 4. Starting Reading
- 4-1 Click the Start Scan button to start reading.



Displays information related to the contents of the currently displayed scan.

Displays the scan results of the 1st, 2nd, 3rd, and 4th scan, starting from the left.

Click on the thumbnails to switch the display.

The scanned area is displayed in the real-time window, as shown below.

The stage is read from the left towards the right.



Reading may be finished at any time before the whole scanning area has completed reading. Click the <u>Stop</u> button when you want to finish reading.

Note1:

If you click Stop during reading, the part that has not been read yet will be saved as an image with a data value of 0 (light intensity of 0).

Note2:

When you click the Stop button, reading itself is canceled. You cannot start reading again from the location where reading stopped.

4-2 When you want to change the display parameters of the real-time window, refer to the explanations below and make settings.



oid: The sigmoid tone curve is used to adjust gradations.



: Drag the adjuster.

You may adjust the density of the read image.

Data of a lower light intensity than the line on the left (Low value) will be displayed as a completely white image, and data of a higher light intensity than the line on the right (High value) will be displayed as a completely black image.

Auto Range Scope : If Auto Range Scope is checked, the Image Reader automatically corrects the optimum tone.

Note:

The Image Reader converts data read from samples to images that have an information of 65536 tones, with 0 being the value for white, and 65535 being the value for black.

The tones are indicated by the horizontal axis of the tone curve graph.

	Magnification: 10:1 (10%) 10:1 (10%) • 5:1 (20%) • 2:1 (40%) • addition, after reading, if you click the display area after click • area after click • • •
4-3	Save as : Save the data with a different file name.
	Launch : Launch the registered analyzing software to display the image.

4-4 To read another fluorescent sample continuously, carry out reading by following the above procedures.

Click the Return button to return to the first Reader Settings window.

Do not open the stage door of the FLA-7000 until the stage has completely returned. If it is opened, close it immediately.

		anu	Eddneri
buttons become active, but the	S Return	bu	itton is grayed out
until the stage has completely r	eturned.		

4-5 Finish reading.Before turning off the power of the FLA-7000, shut down the Image Reader software.

Part **4**

Reading Digitized Samples

- 1. Setting the Reading Conditions
- 1-1 Click the Digitizing button.

The Reader Settings window for the Digitizing mode is displayed.



I Top : Return to the main window from the Digitizing mode.

- 1-2 Make these settings before reading digitized samples. Refer to the following explanations of reading conditions when making settings.
 - Image Folder: Specify where to save the file after reading.
 Click the browse button and specify where to save the file after reading.
 File Name: : : Input the name of a file for saving data of the image.
 You may not start reading unless you input a file name.
 - ③ Comment: :

The comment is saved with the image as a file, and can be viewed with the analyzing software. Input it if necessary.



Note:

To use the Digitizing mode, the LD 473 nm laser must be loaded and the Y520 filter must be registered, or the SHG 532 nm laser must be loaded and the O580 filter must be registered. If these conditions are not fulfilled, the Reader Condition Settings window for the Digitizing mode cannot be accessed.

- (5) PMT : You may set the voltage to be applied to the photomultipliers (PMT) as an integral value within the predetermined range.

200	100	50	25 µm
Short	Reading	g Time	Long
Small	Image F	ile Size	Large

The dynamic range that can be detected is bigger with L5 than with L4. If the signals of the sample are in the L4 range, the density gradation is represented more finely if L4 is selected.



Drag and select the scanning area on the FLUOR stage.



- Load Condition... : Use this button when recalling reading conditions saved with Save Condition... .

Note:

When starting the Image Reader, the settings information from the previous session is displayed.

2006. Mar

- 2. Setting a Digitized Sample on the FLUOR Stage
 After setting the digitized sample on the FLUOR stage, place the fluorescent plate for digitizing on top of it.
 For instructions on setting the digitize sample, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.
- 3. Setting the FLUOR Stage
on the FLA-7000Set the FLUOR stage on the FLA-7000.Set the FLUOR stage on the FLA-7000For instructions on setting the FLUOR stage on the FLA-7000, see the
Fluorescent Image Analyzing System FLA-7000 Operation Manual.
- 4. Starting Reading
- 4-1 Click the Start Scan button to start reading.

Image Reader For FLA-7000	
Reader Scanning Condition	
Save File : TEST	
	Low: 0 High: 65535 Curve: Exponential 💌
Magnification : 10:1 (10%) _ 및 (원,) Q,	∽ Auto Range Scope
	Stop
Save as Status : Method : [Digit473] Laser : LD473	
Launch Filter : 🛅 [Y520]	

The scanned area is displayed in the real-time window, as shown below.

 Next Hoder for 116/1980

 Reader Scanning Condition

 Save File : TEST

 Image: State Scanning Condition

 State Scanning Condition

The stage is read from the left towards the right.

Reading may be finished any time before the whole scanning area has completed reading. Click the stop button when you want to finish reading.

Note1:

If you click Stop during reading, the part that has not been read yet will be saved as an image with a data value of 0 (light intensity of 0).

Note2:

When you click the Stop button, reading itself is canceled. You cannot start reading again from the location where reading stopped.

4-2 When you want to change the display parameters of the real-time window, refer to the explanations below and make settings.

Curve :	Exponential 💌	: Select the ty	pe of tone curve from the pull-
🔽 Auto	Exponential	down menu.	
	Linear	Exponential	The exponential tone curve is
	Sigmola		used to adjust gradations
		Linear:	The linear tone curve is used
			to adjust gradations.
		Sigmoid:	The sigmoid tone curve is
			used to adjust gradations.



High : 50000

Low : 20000

: Drag the adjuster.

You may adjust the density of the read image.

Data that is lighter than the line on the left (Low value) will be displayed as a completely white image, and data that is darker than the line on the right (High value) will be displayed as a completely black image.

Auto Range Scope : If Auto Range Scope is checked, the Image Reader automatically corrects the optimum tone.

Note:

The Image Reader converts data read from samples to images that have an information of 65536 tones, with 0 being the value for white, and 65535 being the value for black.

The tones are indicated by the horizontal axis of the tone curve graph.



Launch

: You may change the display area by selecting a magnification ratio from the pull-down menu. In addition, after reading, if you click the display area after click 🔍 🔍 the buttons, the clicked area can be enlarged or reduced.

4-3 Save as .. : Save the data with a different file name.

: Launch the registered analyzing software to display the image.

- 4-4 To read another digitized sample continuously, carry out reading by following the above procedures. Click the S Return button to return to the first Reader Settings window. Do not open the stage door of the FLA-7000 until the stage has completely returned. If it is opened, close it immediately. When Save as .. scanning finishes, the and Launch buttons become active, but the S Return button is grayed out until the stage has completely returned.
- 4-5 Finish reading.
 Before turning off the power of the FLA-7000, shut down the Image Reader software.

2006. Mar

Part 5

Lasers and Filters, Other Settings

1. Registering Laser and **Filter Combinations** (Method)

Note:

Method settings are registered by a serviceman upon installation. Under normal circumstances, it is not necessary to register these settings.

You may register, change, and delete laser and filter combinations.

1-1 Click the 🖳 Method button on the main window. The following dialog box appears.



Click the Add button. 1-2 The following dialog box appears.

Method	?×	Name:
		Input a
Name : no name		to be re
Laser: LD785	-	Laser
Filter : 1520DF20]	•	Select
	<u>O</u> K <u>C</u> ancel	may al
		I not act

Input a name for the combination to be registered.

Laser :

Select the type of laser. (You may also select lasers that are not actually loaded.)

Filter :

Select the type of filter. (You may also select lasers that are not actually loaded.)

However, combinations of lasers and filters that are not loaded cannot be selected in the Reader Condition screen.

1-3 Input a name for the combination, select the type of laser and filter, and click the <u>uk</u> button. The laser and filter combination is registered .

Note:	Method Setting
Click the Delete and Edit buttons to delete or change a registered Method. Methods initially registered as default cannot be deleted or edited. The Methods marked with [**] are the Methods that are initially set as default.	Mode Settings Method List Method Vame Laser Digitalization Loaras Digitalization SNG522 SYRR0 Ruby LOras Digitalization Loaras Digitalization Loaras Digitalization Loaras Digitalization Loaras Manedoylor LOras Digitalization Loaras Digitalization Loaras Digitalization Loaras Digitalization Loaras Cancel O K

1-4 Click the OK button.

- 2. Filter Module Settings (Filter Module)
 - 2a. Registering Filters in the Image Reader

Note:

Filter Module settings are registered by a serviceman upon installation. Under normal circumstances, it is not necessary to register these settings.

After exchanging the filters of the FLA-7000, you must register the exchanged filters in the Image Reader.

Note:

If filters are not registered in the Image Reader, they are not displayed in the Image Reader window, even if they are physically set in the FLA-7000.

The following explains the method for registering filters in the Image Reader when the [605DF40] is set in filter module position No.2 (second from the back).



i Use this button to change the name or color of the displayed icon for the registered filter. Default filters ([**]) cannot be edited.

- Delete : Use this button to delete a registered filter. Default filters ([**]) cannot be deleted.
- 2a-2 Exchange the filter of the FLA-7000.

For instructions on exchanging the filter of the FLA-7000, see the Fluorescent Image Analyzing System FLA-7000 Operation Manual.

2a-3 Click filter position No. 2.

A red frame appears around the selected filter position.



2a-4 Select [605DF40] from the Filter List. The selected item is highlighted in blue.

Filter Setting		
Filter Module Settings Filter List Name (\$750F20) (\$750F20)	Fifter Module : Test Module 1: [P] 2: [YS20] 3: [S20DF20] 4: [S75DF20]	Load Filter Mudule
Cancel		ОК

2a-5 Click the button.

Filter Setting	
Filter Module Settings	
Fiter List Name / (5750F20) (5750F20) (5750F20) (7570F20) (7570F20) (7570F20) (7570F20) (7530	Filter Module : Test Module •3 BSSDF4H •2 DP •1 Filter Module : Test Module : DP : 0 00507401 : : 0 :
Cancel	ОК

It is also possible to drag-and-drop the selected filter to filter position No. 2.

Filter position No. 2 changes to [605DF40].





2b. Saving the Four Types of Filter Combinations You may save four types of filter combinations that are currently displayed.

Note:

Exchanging the module and saving/recalling filter combinations can be managed more easily if each user has their own filter module.

- 2b-1 Click the Filter Module button on the main window.
- 2b-2 Click the Save Filter Mudule... button. The following dialog box appears.

Save Filter Module	?×
Filter Module Name :	
Test Module	
	·····
	<u>OK</u> <u>C</u> ancel

Filter Module Name:

Input a name for the filter combination.

2b-3 After inputting a name for the filter combination, click the button.

The filter combation is saved.

2b-4 Click the OK button.

- 2c. **Recalling Four Types of Filter Combinations**
- You may recall four types of filter combinations that are currently registered.
- Click the Filter Module button on the main window. 2c-1
- 2c-2 Click the Load Filter Mudule... button. The following dialog box is displayed.

Name		
123	 	
Test Module	 	
Test Module 1		
Test Module 2		
Test Module 3		

- Select the name of the filter combination you want to recall, and 2c-3 click the <u>o</u>^K button. The filter combination is recalled.
- ОК Click the button. 2c-4
- 2d. Registering a New Filter Name

You may register a filter name.

- 2d-1 Click the Filter Module button on the main window.
- 2d-2 Click the Add button. The following dialog box appears.

Silter	Name:
	Input a
Name : no_name	
lcon : 1 -	Icon:
	Select
<u>O</u> K	be disp

put a name for the filter.

elect the color of the filter icon to e displayed in the software.

<u>0</u>K

<u>C</u>ancel

2d-3 Select the filter name and icon color you want to register, and click the 📴 button.

The filter is newly registered.

			Not T fil M	<i>e:</i> o delete or edit the register filter name, select the ter name and click the <u>Delete</u> or <u>Edit</u> button. lethods initially registered as default cannot be eleted or edited.
			2d-4	Click the OK button.
 3. Other Settings (Preference) 3. Other Settings (Preference) 3. Settings for Scanning (Scan Settings) 3. Settings for Scanning (Scan Settings) 3. Click the <u>Preferences</u> button on the main window. 3. Click the <u>Son Settines</u> tab. The following dialog box appears. 		are settings that are carried out when specifying settings for reading es. ding on the reading mode, there may be some functions that cannot d.		
		3a-1	Click the Preferences button on the main window.	
		3a-2	Click the Scan Settings tab.	
		The following dialog box appears.		
				Fluorescence Stage

3a-3 For each of the following items, refer to the explanations below and click the radio button.



When detecting a very small sample amount, the Log format is more effective, because the low-density areas are converted to the gray scale more finely in the Log format. Log format is also recommended for reading gels stained with CBB or silver in the Digitizing mode.

The image will have a low background and clear differences in density. If there are large sample amounts in the IP, Fluorescence, or Digitizing modes, the Root format is more effective.

Note:

The following software versions are required for quantitative analysis of images read in the Root format. Root format files cannot be opened properly with older software versions. For Windows: Science Lab 2003 (Multi Gauge Ver.2.1, Colony Ver.1.1, L-Process Ver.2.1) or later versions For Macintosh: Science Lab 2003 (Image Gauge Ver.4.2, L-Process Ver.2.1) or later versions



Uses specific image shading correction data that was adjusted in accordance with each laser in the Fluorescence mode.

Correction settings will be added in the Fluorescence mode. You will be able to select from optional image shading correction data.

Note:

C On

Registration of the correction settings list is performed by a serviceman. Please contact the dealer where you purchased the FLA-7000, or contact Fuji Photo Film.



Does not use ND filter for adjusting light intensity.

Uses ND filter for adjusting light intensity.

Scan Mode Standa	If select Quick, reading time would become shorter. However, noise would stand out during reading. The reading time varies depending on the setting image size and the scanning area. In case of reading the whole area of FLUOR stage, the reading time are as follows:
	< Standard Mode > 200um : 210sec, 100um : 210sec, 50um : 330sec, 25um : 450sec
	<quick mode=""> 200um : 150sec, 100um : 150sec, 50um : 210sec, 25um : 330sec</quick>
3a-4	A function for supporting new stages in the future.

Ekuwananan Stara	
------------------	--

3a-5 Click the <u>⊆</u> button.

3b. Selecting the File Format for Saving and the Analyzing Software to Launch (Image File Settings) These are settings for saving images. Depending on the scanning mode, there may be some functions that cannot be used.

- 3b-1 Click the References... button on the main window.
- 3b-2 Click the Trace File Settines | tab. The following dialog box appears.

Preferences	
Scan Settings Image File Settings File Format ⓒ img/inf ⓒ img/inf + TIFF Launch Application	
C:¥Program Files¥FUJIFILM¥Science Lab 2005¥MultiGauge.exe Select.	
QKQancel	1

3b-3 Refer to the explanations below and click a radio button.

File Format	The standard file format of the Fuji Film BAS/FLA series. Each file is saved as a combination of an information file (inf file) and a raster file (img file). The Fuji Film Science Lab or Array Gauge may be used for analysis.
	In combination with an img file and inf file, a read image can also be saved in TIFF file format. For TIFF files, image data type is always set to Linear format.

3b-4 Click the <u>select</u> button, and select the specified analyzing software.

Launch Application	
C#Program Files#FUJIFILM#Science Lab 2005#MultiGauge.exe	•
	Select

2006. Mar

Part 6

Installing and Uninstalling the Software

- 1. Installation (For Windows®)
 - 1a. Installation of FUJI USB Control Driver

It installs in the following sequence.

- 1a. Installation of FUJI USB Control driver
- 1b. Installation of FUJI USB Function driver
- 1c. Installation of Image Reader FLA-7000 software

Note:

The computer and FLA-7000 should not be connected with a USB cable during the operation.

- 1a-1 Open the control panel and click "Printers and Other Hardware".
- 1a-2 Click the "Add Hardware".



1a-3 Click the "Next" button.



1a-4 Select "Yes, I have already connected the hardware" and click the "Next" button.



1a-5 Select "Add a new hardware device" and click the "Next" button.



1a-6 Select "Install the hardware that manually select from a list [Advanced]" and click the "Next" button.

Add Hardware Wizard
The wizard can help you install other hardware
The wizard can search for other hardware and automatically install it for you. Or, if you know exactly which hardware model you want to install, you can select it from a list.
What do you want the wizard to do? Search for and install the hardware automatically (Recommended) Install the hardware that I manually select from a list (Advanced)
< <u>₿ack</u> <u>Next</u> Cancel

1a-7 Select "Show All Devices" and click the "Next" button.



1a-8 Click "Have Disk..." button.

Add Hardware Wizard			
Select the device driver you want to install for this hardware.			
Select the manufacturer and model of your hardware device and then click Next. If you have a disk that contains the driver you want to install, click Have Disk.			
Mundahara 🔊 Mada			
(Standard Infrared Port) Serial Cable using IrDA Protocol			
(Standard modern rypes) (Standard port types) (Standard system devices)			
This driver is digitally signed.			
Tearing why unversignmut is induced.			
< <u>B</u> ack Next> Cancel			

1a-9 Click "Browse..." button.



1a-10 Set the folder in which a file is saved to FLA-7000 CD-ROM.



1a-11 Open the FUJIUSB control folder.



1a-12 Select a DevMng.inf file and click the "Open" button.

Locate File			?×
Look jn: 🛅	FUJIUSB Control	Y (3 🗇	⊵•
DevMng.in	f]		
File <u>n</u> ame:	DevMng.inf	*	Open
Files of type:	Setup Information (*.inf)	~	Cancel

1a-13 Click the "OK" button.



1a-14 Click the "Next" button.

Add Hardware Wizard	
Select the device driver you want to in	stall for this hardware.
Select the manufacturer and model of yu have a disk that contains the driver you	our hardware device and then click Next. If you want to install, click Have Disk.
Model Fujifilm USB Device Managment Driver	
This driver is not digitally signed! Iell me why driver signing is important	Have Disk
	< <u>B</u> ack <u>N</u> ext> Cancel

1a-15 Click the "Next" button.



1a-16 Click the "Continue Anyway" button.



1a-17 Click the "Finish" button.



- 1b. Installation of FUJI USB Function Driver
- 1b-1 Connect the USB cable and turn ON the power switch of the FLA-7000.

Note:

Perform this operation, or the personal computer may be reset.

1b-2 Click the "Next" button.



1b-3 Click the "Finish" button.



- 1c. Installation of Image Reader FLA-7000 Software
- 1c-1 Insert the installation CD-ROM of FLA-7000.
- 1c-2 Click the setup icon.
- 1c-3 Click the "Next" button.



1c-4 Click the "Install" button.



1c-5 Click the "Finish" button.



Cancel

- 2. Uninstallation (For Windows[®])
- 2-1 Open the control panel and select "Add or Remove program".



2-2 Select Image Reader FLA-7000 and click the "Remove" button to be deleted.

A progress bar is displayed and uninstallation is started.

Note:

Connection files created after installation of the reader software are deleted. Correction files such as shading data is required to read by FLA-7000. They arelocated in Data folder of FLA-7000 folder. The files you wish to keep must be saved outside the FLA-7000 folder.

3. Installation To install Image Reader FLA-7000 software, follow the procedure described below. (For Macintosh[™]) below.

Installation of FLA-7000

< For Mac OS10.4 >

- 3-1 Double-click the "FLA-7000 Install CD" icon to open it.
- 3-2 Click the "FLA-7000.pkg" icon.
- 3-3 Click the "Continue" button on the "Introduction" screen.
- 3-4 Click the "Continue" button on the "Installation Destination" screen.
- 3-5 Click the "Upgrade" button on the "Installation Type" screen.
- 3-6 Provide the name and password of the administrator.
- 3-7 Click the "close" button on the "Finish Up" screen.
- 4. Uninstallation
(For Macintosh™)Move the FLA-7000 folder contained in the Application folder on the Machintosh
HD to your Trash bin.

Note:

Connection files created after installation of the reader software are deleted. Correction files such as shading data is required to read by FLA-7000. They arelocated in Data folder of FLA-7000 folder. The files you wish to keep must be saved outside the FLA-7000 folder.

Part 7

Troubleshooting

1. Errors

An error is the condition in which all of the reading modes of the FLA-7000 cannot be used.

When a warning occurs, the indicator lamps make a sound ("Be-beep, be-beep, be-beep").

An error message dialog box is displayed on the Image Reader screen displayed on the computer.



Please contact Fuji Photo Film with the 4-digit Error Code and the 8-digit number inside the parentheses.

2. Warnings

The following are examples of warnings that are displayed as a message in the Image Reader window. If a displayed message includes instructions, please follow them.

Message	Meaning and Countermeasure
Filter module has not been setup.	Meaning:
Push the filter module button and	• A filter module has not been set.
set up the filter module.	• Press the filter module button, and set a filter module.
	Countermeasure:
	Confirm that a filter module is set.
The combination of the lasers and filters could be inappropriate. Check the lasers and filters.	Meaning: There is a possibility that the laser and filter combination is inappropriate. Please check the laser and filter. Countermeasure: • Confirm the Methods registered. It is possible that the laser and filter combination is inappropriate. • Confirm that the filter specified in the Image Reader is actually set.

Laser error occurred. Use other lasers.	 Meaning: A laser error has occurred. Use a different laser. Countermeasure: An error has occurred with the laser. If an error occurs even after restarting the instrument, call a serviceman. It is possible to scan using a different laser.
Failed to retrive picture data from PC. Check the PC setting.	 Meaning There was a failure in reading the image data from the computer. Check the computer environment. Countermeasure: It is possible that the computer's processing ability has degraded. Check the computer once, when there are no devices connected to it and no software started.
PMT error occurred. Use low sensitivity for the setting.	 Meaning: An overexposure error has occurred. Scan with low sensitivity settings. Countermeasure: Overexposure has been detected outside of the scanning area. If an error occurs even after restarting the instrument, call a serviceman.
The file name is already used. Close the SampleSetDoor or FilterChangeDoor.	Meaning: Close the sample set door or the filter change door. Countermeasure: Close the sample set door or the filter change door.

Note:

If an error message is displayed, a serviceman should take the countermeasures to resolve the trouble. Please contact the dealer where you purchased the FLA-7000, or contact Fuji Photo Film.

Fuji Photo Film Co., Ltd. 2-26-30 Nishi Azabu, Minato-ku Tokyo, Japan 106-8620 TEL: +81-3-3406-2201 FAX: +81-3-3406-2158 e-mail: sghelp@fujifilm.co.jp

For technical questions, send an e-mail to the following address. e-mail: sginfo@fujifilm.co.jp.